'PRODUCT NAME DCR BODKINS STRAIGHT AND ANGLED
'PRODUCT INFORMATION AND INSTRUCTIONS FOR USE

Symbols Used

LOT				STERILE R			EC REP
GB- Lot Number ES - Número de lote FR- N° de lot IT - Numero di lotto	GB- Use Until Date ES - Usar antes de la fecha FR- Date limite d'utilisation IT - Utilizzare entro	GB- Do Not Re-use ES - No reutilizar FR- Ne pas réutiliser IT - Non riutilizzare	GB- Do not use if product is opened or damaged. ES - No utilizar si el producto está abierto o dañado FR- Ne pas utiliser si l'emballage est ouvert ou endommagé. IT - Non usare se la confezione è aperta o danneggiata	GB- Sterilised by Irradiation ES- Esterilizado por irradiación FR- Stérilisé par irradiation. IT- Sterilizzato a raggi gamma.	GB - Manufacturer ES - Fabricante FR - Fabricant IT - Fabbricante	GB - Date of Manufacture ES - Fecha de fabricación FR - Date de fabrication IT - Data di fabbricazione	GB - EU Authorised Representative ES - Representante autorizado de la UE FR - Représentant UE autorisé IT - Mandatario nell'Unione europea
REF					MD	UDI	
GB- Catalogue Number ES - Número de catálogo FR- Numéro du catalogue IT - Riferimento di catalogo	GB- Caution ES - Advertencia FR- Avertissement IT - Attenzione	COR_016 GB - Consult Electronic Instructions for Use ES - Consultar instrucciones de empleo electrónicas FR - Consulter le mode d'emploi électronique IT - Consultare le Istruzioni per l'uso elettroniche	R _{only} GB - Caution: United States Federal law restricts medical devices to sale by or on the order of a physician ES - Advertencia: las leyes federales de Estados Unidos limitan la venta y el uso a un facultativo (o por prescripción de este) FR- Avertissement : La loi fédérale américaine limite la vente et l'utilisation aux médecins ou à la demande d'un médecin IT - Attenzione: la legge federale degli Stati Uniti d'America limita la vendita e l'uso ai medici o su prescrizione del medico	GB- Do not resterilise ES - No reesterilizar FR- Ne pas restériliser IT - Non risterilizzare	GB - Medical Device ES - Dispositivo medico FR - Dispositif médical IT - Dispositivo medico	GB - Unique Device Identifier ES - Identificador único del dispositivo FR - Identifiant unique du dispositif IT - Identificatore univoco del dispositivo	GB - Single sterile barrier system with protective packaging outside ES- Sistema de barrera estéril individual con embalaje protector en el exterior FR - Système de barrière stérile unique avec emballage protecteur externe IT - Sistema di barriera singola sterile con confezionamento


Fig. 1
GB - This Instructions for use relates to the following products:
ES - Estas instrucciones de uso hacen referencia a los siguientes productos.
FR - Mode d'emploi des produits suivants :
IT - Le seguenti istruzioni per l'uso si riferiscono ai seguenti prodotti:


Fig.2


Fig.3

GB	DCR Bodkin Angled	DCR Bodkin Straight
ES	DCR Bodkin acodado	DCR Bodkin recto
FR	Sonde de Bodkin coudée pour DCR	Sonde de Bodkin droite pour DCR
IT	Sonda DCR Bodkin angolata	Sonda DCR Bodkin dritta

GB - INSTRUCTIONS FOR USE

Description

NETWORK MEDICAL DCR Bodkins consist of two malleable stainless steel needles attached by a length of silicone rod. Devices are supplied sterile for SINGLE USE ONLY.

Intended Use

DCR Bodkins are intended for probing and /or silicone intubation during a dacryocystorhinostomy (DCR) procedure. Dacryocystorhinostomy (DCR) surgery is a procedure that aims to eliminate fluid and mucus retention within the lacrimal sac, and to increase tear drainage for relief of epiphora (water running down the face). A DCR procedure involves removal of bone adjacent to the nasolacrimal sac and incorporating the lacrimal sac with the lateral nasal mucosa in order to bypass the nasolacrimal duct obstruction. This allows tears to drain directly into the nasal cavity from the canaliculi via a new low-resistance pathway.

CAUTIONS

- This device is supplied STERILE and ready to use.
- This device is for SINGLE USE ONLY. Do NOT re-sterilise or re-use.
- Do not use if the packaging has been opened or damaged.
- This device is only to be used by a suitably trained or qualified person under the facility approved and recommended conditions.
- Ease of passage of the probes through the Lachrymal System will vary widely from patient to patient with the placement of the tubes in some patients being quite difficult due to narrow openings in the lachrymal system, especially at the lower end of the bony canal. Separation of the silicone tubing from the probe may occur when intubating patients with severe blockage or a pronounced bony prominence.
- CAUTION the silicone rod should be implanted for no more than 29 days.**
- Silicone tubing could be difficult detect or to remove. Examination under anaesthesia, including nasal endoscopy may be warranted when accidental tube extrusion cannot be reliably confirmed.
- For safe disposal, the needle should be considered as 'sharp' and must be immediately placed in a secure puncture-resistant bin (conforming to BS7320) which is suitable for incineration. All the components should be disposed of in accordance with hospital approved procedures.
- The product comes into contact with bodily fluids, which can be contaminated. Care should be taken in the handling and disposal of the silicone tubing after use to prevent contamination.
- Do not use for any purposes other than lacrimal surgery.

MRI Compatibility

All tubing made from Silicone, which is MRI compatible.

Incident Reporting

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established

Instructions and Surgical Procedure

- Figure 1 illustrates the anatomy of the lacrimal system.

**'PRODUCT NAME DCR BODKINS STRAIGHT AND ANGLED
'PRODUCT INFORMATION AND INSTRUCTIONS FOR USE**

- Each punctum should be dilated.
- During the DCR procedure, a small incision is made either in the skin or inside the nose.
- A small piece of bone is removed.
- The bodkin probes, with silicone attached, are inserted through both upper and lower puncta, fed through the canaliculi and orientated down through the lachrymal system into the nose through the ostium created (fig 2), approximately 35mm - 40mm.
- The silicone rods are tied (a suture may be used for this) and left long in the nasal vestibule to facilitate removal later (fig 3).
- The lacrimal sac and the nasal mucosa are brought together and the incision is closed.
- The silicone rod is left in position as a stent, which can remain for up to 29 days.
- Removal of the rod is accomplished by locating the rod between the puncta and pulling the silicone upwards. The rod can then be cut and removed.
- The silicone rod must be removed through the nose.

Possible Adverse Effects

- The rod can be implanted for up to 29 days, however the surgeon should regularly monitor the patient for suitability. The implant should be removed or replaced if damaged when deemed appropriate by the surgeon.
- Complications which may occur include, but are not limited to: granuloma formation, cellulitis, infection, epistaxis, punctal or canaliclar erosion, corneal erosion or abrasion, asymptomatic cheesewiring or masupialization, nasal discharge, periorbital edema and periorbital ecchymosis.
- It is important to note that Canaliculi or puncta may become lacerated or split if the tubes are too tight. Care must be taken in children to allow for patient growth. The surgeon should regularly monitor the patient for any signs of damage around the implant for punctal erosion, elongation or any sign of tube migration through the puncta.
- Early post-surgical displacement of the implant may occur particularly in children or when forcibly sneezing, potentially leading to lid notching or development of an ectropion due to poor healing.
- Patients requiring Continuous Airway Pressure (CPAP) therapy may experience irritation following the procedure due to air regurgitation.
- Particular attention should be taken during and prior to knotting/suturing to avoid over extension of tubing.
- There is a risk of synechia formation between the lateral nasal wall and the middle turbinate mucosa.
- Tubing is not X-Ray or CT Scan detectable.
- The probes of product contain nickel which may cause allergic reaction in patients with sensitivity to nickel. Risk assess the use of the devices in relation to the medical benefit of the procedure and take necessary precautions with patients with known sensitivity to nickel.
- This product contains silicone which may cause allergic reaction in patients with sensitivity to silicone. Risk assess the use of the devices in relation to the medical benefit of the procedure and take necessary precautions with patients with known sensitivity to silicone.
- There is a potential risk of synechia formation between the lateral nasal wall and the middle turbinate mucosa
- There is a potential risk of cross infection if devices are not disposed of after use following hospital approved procedures.
- There is a risk of periorbital edema and/or periorbital ecchymosis.

Sterilisation

The NETWORK MEDICAL DCR Bodkins are supplied sterile by Gamma Irradiation and are designated for SINGLE USE ONLY.

HAZARDS ASSOCIATED WITH THE RE-USE OF SINGLE USE ONLY DEVICES

1. Single use devices have not been validated for re-use.
If you re-use a device you may be held Legally Liable for the safe performance.
2. Cross-contamination and infection risks to patients. Including transmission of:
 - CJD & Variant CJD.
 - Prion Diseases.
 - Bacterial Endotoxins.
 - Hepatitis B & Hepatitis C.
 - Risks posed by HIV and AIDS
3. Device failure through material fatigue or degradation caused by initial use and design:
 - Plastics: Can be weakened, warped or become brittle.
 - Metals: Can be damaged or subjected to rusting. Blade edges will be damaged or blunt. Needles may be bent from the original shape.
 - Adhesives: Can degrade if re-sterilised causing parts to dislodge.
4. Patient injury from device failure and/or chemical burns from residue of decontamination agents absorbed into the materials.

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