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PRODUCT NAME: EPISTAXIS CATHETER
PRODUCT INFORMATION AND INSTRUCTIONS FOR USE

NL	Epistaxis-katheter groot	Epistaxis-katheter klein
PT	Cateter para Epistaxe Grande	Cateter para Epistaxe Pequeno
SV	Epistaxiskateter stor	Epistaxiskateter liten

GB - INSTRUCTIONS FOR USE

Description

The Network Epistaxis Catheter is an all silicone catheter with two independently inflatable cuffs. The airway permits nasal breathing and access for suctioning.

The contents of this package are sterile unless the package is opened or damaged.

Intended Use

Epistaxis catheters are used to stem serious posterior and/or anterior nasal bleeding (epistaxis), which cannot be stemmed by other methods. They establish tamponade by traction of the balloon.

CAUTIONS

- This device is supplied STERILE and ready to use.
- This device is for SINGLE USE ONLY. Do NOT re-sterilise or re-use.
- Do not use if the packaging has been opened or damaged.
- This device is intended for use by trained medical persons possessing the requisite skill and experience to use the device in accordance with the prevailing standards of medical practice and in conjunction with the instructions for this device.
- The product comes into contact with bodily fluids, which can be contaminated. Care should be taken in the handling and disposal of the device after use to prevent contamination.
- CAUTION: US Federal law restricts this device to sale by or on the order of a physician
- Do not use petroleum based lubricants with this device as it will seriously affect the performance of the balloons potentially causing them to burst.
- Sterile saline is not to be used to inflate the device as the salts used may crystallise and cause a blockage.
- This product contains silicone which may cause allergic reaction in patients with sensitivity to silicone. Risk assess the use of the devices in relation to the medical benefit of the procedure and take necessary precautions with patients with known sensitivity to silicone.
- This device has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the device in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.
- This device is compatible with a syringe which meets the requirements of BS EN 20594-1:1994 (ISO 594-1:1986).
- Care must be taken not to over-inflate balloons.

Incident Reporting

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established

INSTRUCTIONS

- Pre-use checks:
 - Attach a syringe to the pilot port marked 10 and test-inflate the distal cuff with 5 to 10 ml of air. Deflate cuff by depressing the valve using a syringe. Evacuate the cuff completely and when balloons appear collapsed gently withdraw the catheter.
 - Attach syringe to the pilot port marked 30 and test-inflate the proximal cuff with 15 to 25 ml of air. Deflate cuff by depressing the valve using a syringe. Evacuate the cuff completely and when balloons appear collapsed gently withdraw the catheter.
 - If either cuff does not inflate or fully deflate, use a new tube. Save the old tube and call Network Customer Service.
- Lubricate a catheter using a water-soluble lubricant and insert the catheter so that the bevelled tip is in the nasopharynx.
- Attach a syringe to the pilot port marked 10 and inflate the distal cuff with air or sterile fluid until the desired seal is obtained. Volumes of 4 to 8 ml are usually adequate. **Do not over-inflate or exceed 10 ml.** Remove the syringe and gently pull the catheter anteriorly to seat it.
- Attach a syringe to the pilot port marked 30 and inflate the proximal cuff with air or sterile fluid until the desired seal is obtained. Volumes of 10 to 25 ml are usually adequate. **Do not over-inflate or exceed 30 ml.** Remove the syringe.
- Monitor the cuff pressures frequently and adjust as needed.
- A small well lubricated suction catheter may be inserted through the airway to remove secretions and/or clots.

Removal

Evacuate each cuff completely by depressing the valve using the syringe and when balloons appear collapsed gently withdraw the catheter.

Sterilisation

- The device is a SINGLE USE ONLY device supplied sterile and ready for use. Sterilisation is by Ethylene Oxide (EO).

HAZARDS ASSOCIATED WITH THE RE-USE OF SINGLE USE ONLY DEVICES

- Single use devices have not been validated for re-use.
If you re-use a device you may be held **Legally Liable** for the safe performance.
- Cross-contamination and infection risks to patients. Including transmission of:
 - CJD & Variant CJD.
 - Prion Diseases.
 - Bacterial Endotoxins.
 - Hepatitis B & Hepatitis C.
 - Risks posed by HIV and AIDS
- Device failure through material fatigue or degradation caused by initial use and design:
 - Other materials:** May degrade, becoming unacceptable when compared to original manufacturing criteria.
- Patient injury from device failure and/or chemical burns from residue of decontamination agents absorbed into the material

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PRODUCT NAME: EPISTAXIS CATHETER

PRODUCT INFORMATION AND INSTRUCTIONS FOR USE

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