REF EP-3000

Middle Ear Inflation Device
Appareil d'expansion de l'oreille moyenne
Inflationsgerät für das Mittelohr
Dispositivo de inflado del oído medio
Dispositivo per il gonfiaggio dell'orecchio medio
Instrument voor middenoorinsufflatie
Oppblåsningsutstyr for mellomøre

Manufactured by:

815 Vikings Parkway, Suite 100
St. Paul, MN 55121 | USA

P: 1-888-229-2875 | 651-789-3939
F: 1-888-229-1941 | 651-789-3979

www.summitmedicalusa.com

Instructions for Use
Mode d'emploi
Gebrauchsanweisung
Instrucciones de uso
Istruzioni per l'uso
Gebruiksaanwijzing
Bruksanvisning

INDICATIONS FOR USE

The EarPopper is indicated for use in the treatment of negative middle ear pressure.

Negative middle ear pressure can lead to fluid accumulation in the middle ear, impaired hearing and hearing loss. The EarPopper provides a method for ventilating the middle ear by introducing air through the nose and Eustachian tube into the middle ear. Equalizing middle ear pressure can prevent the accumulation of fluid and prevent hearing loss.

PRECAUTIONS

- Do not force the instrument into the nasal cavity.
- To ensure compliance, the physician should demonstrate the proper use of the device when first recommending the EarPopper.
- Do not use the EarPopper without a clean, intact nosepiece attached.

CONTRAINDICATIONS

- Patients who will be using the device at home should be instructed not to use the EarPopper in the presence of severe ear infection, upper respiratory infection, heavy nasal congestion, or cold symptoms without first consulting their physician.

ADVERSE REACTIONS

- Patients may feel a pressure change and/or “pop” in one or both ears after use of the EarPopper. The ears may feel temporarily “clogged” after inflation. This sensation is an indication that the procedure is working correctly and it should decrease within a few minutes.

PARTS OF THE EARPOPPER

See instructional video online at YouTube website:

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OPERATING INSTRUCTIONS (Continued):

2. Hold the nosepiece portion of the EarPopper firmly against one nostril. Press the other nostril closed using the index finger of the other hand. It is very important to establish a good seal. Press the power button to begin airflow, after 2-3 seconds of airflow, instruct the patient to swallow with mouth closed. See Figure 2.

Figure 2

NOTE:

If the patient does not feel any change in the ear, repeat the procedure making sure the EarPopper remains firmly against the nostril, while the patient swallows with mouth closed. Instructing the patient to hold a small amount of liquid in mouth prior to the procedure will aid in swallowing.

3. Repeat the procedure for the other nostril.
4. Wait five minutes and repeat steps 2-3.
5. To remove nose piece, turn device upside-down over a trash receptacle and press the nosepiece release buttons on both sides of the nosepiece simultaneously. See figure 3.

Figure 3

BATTERY REPLACEMENT:

To check for a dead battery use the following steps:

1. Place finger over hole in the end of the nose piece.
2. Press the Power Switch to activate device.
3. If motor runs intermittently or stops, replace the battery.

Open the battery compartment by removing the screw and sliding out the compartment door. Install a new 9 volt alkaline battery. Replace the battery door and re-install the screw. See Figure 4.

Figure 4

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CARE, CLEANING, & STORAGE

WARNING: TO AVOID CROSS CONTAMINATION REPLACE THE DISPOSABLE NOSEPIECE AFTER EACH PATIENT.

CLEANING: The outside of the device may be cleaned by wiping with a soft, damp cloth. DO NOT AUTOCLAVE, OR USE ABRASIVE CLEANERS, OR CHEMICALS, OR DISINFECTANTS. AND DO NOT SUBMERGE THE EARPOPPER DEVICE IN LIQUID OF ANY KIND.

STORAGE: Store the EarPopper in a clean dry location at room temperature with a clean nosepiece in place to protect it from damage and dirt.

Remove battery for prolonged storage (6 months or more).

DEVICE DISPOSAL

The EarPopper contains no lead, cadmium, mercury or bromide flame retardants.

All plastic pieces (body top/bottom and misc. internal plastic parts) can be recycled as "general plastic".

Internal wires, motor and brass can be recycled as "general metals".

Battery - dispose of according to local regulations.

WARRANTY

This device is warranted to be free from manufacturing defects for a period of one year from date of purchase under normal use. This excludes damage to the product resulting from accident or misuse.

CUSTOMER SERVICE

For questions concerning the operation of the EarPopper, or to request accessories and replacement parts, please visit www.earpopper.com or call or write using the contact information on the back of this manual.

TECHNICAL INFORMATION

Product Safety Compliance Testing: (IEC 60601-1)

IEC 60601-1-1

IEC 60601-1-2

EN61000-4-2 Electrostatic Discharge	PASS	8 KV air 6 KV contact
EN61000-4-3 Radiated Immunity	PASS	10 V/M
EN61000-4-4 Electrical Fast Transients	N/A	—
EN61000-4-5 Surge Immunity	N/A	—
EN61000-4-6 Conducted Immunity	N/A	—
EN61000-4-8 Power Frequency Magnetic Field	PASS	1 amp/M
EN61000-4-11 Voltage Dips, Interrupts, and Fluctuations	N/A	—

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