
'PRODUCT NAME CLEAR-VUE™ RECIPIENT VACUUM TREPHINE FOR PK
'PRODUCT INFORMATION AND INSTRUCTIONS FOR USE

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DIAGRAM 1

PRODUCT NAME CLEAR-VUE™ RECIPIENT VACUUM TREPHINE FOR PK

PRODUCT INFORMATION AND INSTRUCTIONS FOR USE

GB - This Instructions for use relates to the following products:
 DE - Diese Gebrauchsanweisung bezieht sich auf die folgenden Produkte:
 ES - Estas instrucciones de uso hacen referencia a los siguientes productos:
 FR - Mode d'emploi des produits suivants :
 IT - Le seguenti istruzioni per l'uso si riferiscono ai seguenti prodotti:
 PT - Estas instruções referem-se aos seguintes produtos:
 RU - Настоящая инструкция по эксплуатации содержит информацию о следующих продуктах:
 TR - Bu kullanım kılavuzu aşağıda belirtilen ürünler ile ilgilidir:
 ZH - 本使用说明涉及以下产品:

GB	Clear-Vue™ Recipient Vacuum Trephine for PK
DE	Clear-Vue™ Empfängerhornhaut-Vakuum-Trepan für PK
ES	Clear-Vue™ Recipient Vacuum Trephine para QP
FR	Trépan à vide receveur Clear-Vue™ pour PK
IT	Trapano corneale a suzione ricevente per PK Clear-Vue™
PT	Trépano de Vácuo para Recetor Clear-Vue™ para PK (penetração ceratoplástica)
RU	Вакуумный трепан для реципиента Clear-Vue™ для проникающей кератопластики
TR	PK için Clear-Vue™ Vakumlu Alıcı Trefini
ZH	用于 PK 的 Clear-Vue™ 受体真空环钻

GB – INSTRUCTIONS FOR USE

Description

The Clear-Vue™ Recipient Vacuum Trephine for PK has an ultra-sharp circular blade held inside an outer suction ring and attached to a spring-loaded syringe. Evacuating the ring when the trephine is placed in contact with the cornea, creates a vacuum which holds the trephine firmly in place on the cornea. Crosshairs aid centration over the optical axis.

Intended Use

The CORONET range of corneal trephines are devices used in Corneal Graft Surgery and can be used for full penetrating keratoplasty of the recipient cornea or for the preparation of the donor button during corneal transplant procedures. The ultra-sharp, thin profile cutting edge minimises damage to the endothelial cells of the cornea and reduces undercut and tissue distortion.

The Clear-Vue Recipient Vacuum Trephine for PK is a single use trephine designed for full penetrating keratoplasty.

CAUTION

- The device is for SINGLE USE ONLY. Do NOT re-sterilise or re-use.
- This device is intended for use by trained medical persons possessing the requisite skill and experience to use the device in accordance with the prevailing standards of medical practice and in conjunction with the instructions for this device.
- This device is supplied sterile and ready for use. Do not use if the pouch has been opened or damaged.
- The Coronet trephine is much sharper than a traditional trephine. Caution should be exercised at all times when handling and using the trephine.
- The Coronet trephine is supplied with the blade in perfect condition. Should the device be accidentally dropped onto a hard surface, check the blade integrity under the microscope before use and discard the blade if damaged.
- The product comes into contact with bodily fluids, which can be contaminated. Care should be taken in the handling and disposal of the device after use to prevent contamination.
- CAUTION: US Federal law restricts this device to sale by or on the order of a physician
- This product contains nickel which may cause allergic reaction in patients with sensitivity to nickel. Risk assess the use of the devices in relation to the medical benefit of the procedure and take necessary precautions with patients with known sensitivity to nickel.
- For safe disposal, the device should be considered as 'sharp' and must be immediately placed in a secure puncture-resistant bin (conforming to BS7320) which is suitable for incineration. All the components should be disposed of in accordance with hospital approved procedures.
- For safety reasons the trephine blade must be aligned with the inner sleeve before placing on the cornea (diagram 1). Retract the blade until it is aligned with the bottom of the inner sleeve i.e. no longer visible.
- Do not retract the blade out of the inner sleeve to prevent damage to the blade edge.

Incident Reporting

- Any serious incident that has occurred in relation to this device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established

Surgical Technique

- The trephine is supplied with one grooved marker pin on the capstan to aid blade control.
- For safety reasons the trephine is supplied with the blade just visible in the inner sleeve (diagram 1).
- Prior to use, retract the blade until it is aligned with the bottom of the inner sleeve i.e. no longer visible.
- Depress the plunger of the suction syringe and fully evacuate the syringe and hold.
- Centre the trephine correctly on the cornea using the crosshair.
- Press down gently on the trephine to seat it on the cornea.
- Release the plunger and wait approximately 5 seconds to allow the vacuum to become fully stable. The piston will travel through the syringe and will stop when vacuum is established.
- Lift the trephine gently to test whether the vacuum has taken and then allow the cornea to return to its normal curvature before rotating the blade.
- Advance the blade to contact with the cornea. The blade travels 63µm (0,063mm) for every quarter turn and 250µm (0,25mm) per complete rotation.
- Continue to advance the blade slowly and proceed with the cut.
- Stop the blade rotation at the first sign of aqueous humour.
- Release the suction by depressing the plunger of the syringe or by cutting the tubing and remove the trephine.
- Complete the procedure in the desired manner.

STERILISATION

- The device is a SINGLE USE ONLY device supplied sterile and ready for use. Sterilisation is by Ethylene Oxide (EO).
- DO NOT RE-STERILISE. Network Medical Products Ltd. assumes no liability for devices that have been re-sterilised by healthcare facilities.

HAZARDS ASSOCIATED WITH THE RE-USE OF SINGLE USE ONLY DEVICES:

- Single use devices have not been validated for re-use. If you re-use a device you may be held **Legally Liable** for the safe performance.
- Cross-contamination and infection risks to patients. Including transmission of:
 - CJD & Variant CJD.
 - Prion Diseases.
 - Bacterial Endotoxins.
 - Hepatitis B & Hepatitis C.
 - Risks posed by HIV and AIDS
- Device failure through material fatigue or degradation caused by initial use and design:
 - Plastics: Can be weakened, warped or become brittle.
 - Metals: Can be damaged or subjected to rusting. Blade edges will be damaged or blunt. Needles may be bent from the original shape.
 - Adhesives: Can degrade if re-sterilised causing parts to dislodge.
 - Other materials: May degrade, becoming unacceptable when compared to original manufacturing criteria.
- Patient injury from device failure and/or chemical burns from residue of decontamination agents absorbed into the materials.

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**'PRODUCT NAME CLEAR-VUE™ RECIPIENT VACUUM TREPHINE FOR PK
'PRODUCT INFORMATION AND INSTRUCTIONS FOR USE**

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