Otologic Ventilation Tubes

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
Paukenröhrchen

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Otologische ventilatiebuisjes

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Instructions for Use

Mode d'emploi

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Bruksanvisning

Otologic Ventilation Tubes

Before using, read the following information:

DESCRIPTION

Otological Ventilation Tubes are small tubular implants available in a variety of biocompatible materials including silicone elastomer, fluoroplastic, polyethylene, and titanium. Numerous designs and sizes are available with single or multiple flanges to satisfy various surgical techniques for insertion and to facilitate communication of the auditory canal with the middle ear.

INTENDED USE

When inserted through a myringotomy, a ventilation tube provides a passageway for movement of air between the auditory canal and the middle ear. The unobstructed passageway may also allow a means of drainage of fluids resulting from acute or chronic otitis media from the middle ear into the auditory canal. In addition, surgical placement of ventilation tubes also provides a means of equalizing air pressure between the outer ear and the middle ear and continued ventilation to prevent fluid accumulation within the middle ear.

INDICATIONS FOR USE

The surgeon must use medical judgment and consider the patient's medical history prior to a decision to surgically insert a ventilation tube. Pathologic conditions for which ventilation tubes are indicated include but are not limited to:

- Eustachian tube dysfunction
- Chronic otitis media with effusion characterized as serous, mucoid, or purulent
- Recurrent acute otitis media which fails to respond satisfactorily to alternative therapies.

CONTRAINDICATIONS

This device is not designed, sold or intended for use except as indicated or prescribed by a physician. Some conditions for which ventilation tubes may not be the therapy of choice include:

- Cases of otitis media which respond favorably to drug therapy alone.

PRECAUTIONS

Avoid fluid contact with ear.

POSSIBLE ADVERSE EFFECTS

- The middle ear may develop subsequent secondary infections from either water or airborne pathogens (Care should be taken to avoid water entering the ear)
- Clogging of the tube, which may require removal and replacement
- Persistent perforation of the tympanum after tube extrusion
- Granulomatous reaction of the tympanum necessitating tube removal
- Premature extrusion of the tube
- Failure of the tube to self-extrude requiring medical intervention for removal
- Allergic reaction may occur in patients with metals sensitivity
- Medialization or dislocation of the tube into the middle ear cavity
- Tympanosclerosis/myringosclerosis
- Localized or diffuse membrane atrophy
- Otalgia, or an ear ache
- Otorrhea, or ear discharge/drainage
- Biofilm formation

WARNINGS

- The surgeon is encouraged to carefully choose a tube size, style and material to meet the needs of each specific clinical situation, while minimizing adverse effects.

MRI INFORMATION

The Stainless Steel Vent Tubes, Titanium Vent Tube and Stainless Steel Wire have been determined to be MR-conditional according to the terminology specified in the American Society for Testing Materials (ASTM) International. Designation: F2503-05.

A patient with this device can be scanned safely immediately after placement under the following conditions:

- Static magnetic field of 3 Tesla or less
- Maximum spatial gradient magnetic field of 720-Gauss/cm or less

Additional information can be obtained by contacting Summit Medical at 888-229-2875.

STERILITY

The metal Otologic Ventilation Tubes are provided STERILE, being terminally sterilized by ethylene oxide (ETO) gas in a controlled environment and confirmed through biological testing. The package should be carefully inspected for punctures, tears and security of the seal or any other evidence of the sealed pouch having been compromised prior to placement of the contents in the sterile field. The product is sterile in an undamaged, unopened package.

The Otologic Ventilation Tubes are single use only and should not be resterilized.

Adhesive product labels are included on the pouch for use on patient and hospital records.

!USA Rx Only

For US Audiences Only

Single-use, one-patient device will degrade if reprocessed. No effective cleaning process has been developed to prevent cross contamination. Contamination of a reprocessed device may lead to injury, illness or death of the patient.

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