


2mm, 3mm, 4mm Suction Tube Range Specification Sheet



Product Overview

The NETWORK ENT® Suction Tube range has been specifically designed to provide precision and control during ENT procedures. Each suction tube is manufactured without bending or pinching to create a continuous radius, significantly reducing the risk of blockage.

- Constructed from high quality Japanese steel for superior thin walled tubing
- Ergonomically balanced and lightweight handle for unrivalled manoeuvrability and control
- Available with and without suction control
- Supplied sterile, single use only, declared 5 year shelf life
- Manufactured in the UK

2mm, 3mm, 4mm Suction Tube Range Specification Sheet

2mm Diameter Suction Tube Size Options

Product Size	Box Quantity	Product Code
Oesophageal/ Laryngeal 2mm x 215mm, 30 degree angle, proximal bend with suction control hole	25	74-2030
Oesophageal/ Laryngeal 2mm x 275mm, 30 degree angle, proximal bend with suction control hole	25	74-2033
Oesophageal/ Laryngeal 2mm x 350mm, 30 degree angle, proximal bend with suction control hole	25	74-2037
Zoellner type, 2mm x 85mm, 10 degree shallow distal bend, without suction control hole	50	74-2011
Zoellner type, 2mm x 85mm, distal bend acute 30 degree angle, without suction control hole	50	74-2012
Zoellner type, 2mm x 85mm, distal bend acute 30 degree angle, without suction control hole	50	74-2012P
Zoellner type, 2mm x 85mm, 10 degree shallow distal bend, with suction control hole	50	74-2013
Zoellner type, 2mm x 85mm, distal bend acute 30 degree angle, with suction control hole	50	74-2014
Zoellner type, 2mm x 85mm, distal bend acute 30 degree angle, with suction control hole	50	74-2014P
Belucci type, 2mm x 85mm, proximal bend with suction control hole	50	74-2015
Zoellner type, 2mm x 100mm, medial bend acute 30 degree angle, without suction control hole	50	74-2016
Zoellner Pen grip type, 2mm x 100mm, medial bend acute 30 degree angle, without control hole	50	74-2016P
Frazier Lampert type, 2mm x 125mm, proximal bend with suction control hole	50	74-2017

3mm Diameter Suction Tube Size Options

2mm, 3mm, 4mm Suction Tube Range Specification Sheet

Product Size	Box Quantity	Product Code
Magill Type 3mm x 110mm, 10 degree angle, shallow distil bend, with suction control hole	50	74-2010
Fraizer Lampert Type 3mm x 125mm, proximal bend, with suction control hole	50	74-2050
Oesophageal/ Laryngeal 3mm x 165mm, 30 degree angle, proximal bend with suction control hole	25	74-2051
Oesophageal/ Laryngeal 3mm x 250mm, 30 degree angle, proximal bend with suction control hole	25	74-2052
Oesophageal/ Laryngeal 3mm x 275mm, 30 degree angle, proximal bend with suction control hole	25	74-2053
Oesophageal/ Laryngeal 3mm x 350mm, 30 degree angle, proximal bend with suction control hole	25	74-2057
Oesophageal/ Laryngeal 3mm x 380mm, 30 degree angle, proximal bend with suction control hole	25	74-2058
Oesophageal/ Laryngeal 3mm x 530mm, 30 degree angle, proximal bend with suction control hole	25	74-2059

4mm Diameter Suction Tube Size Options

Product Size	Box Quantity	Product Code
Oesophageal/ Laryngeal 4mm x 120mm, 30 degree angle, proximal bend with suction control hole	25	74-2077
Oesophageal/ Laryngeal 4mm x 135mm, 30 degree angle, proximal bend with suction control hole	25	74-2081
Oesophageal/ Laryngeal 4mm x 150mm, 30 degree angle, proximal bend with suction control hole	25	74-2079
Oesophageal/ Laryngeal 4mm x 165mm, 30 degree angle, proximal bend with suction control hole	25	74-2080
Oesophageal/ Laryngeal 4mm x 275mm, 30 degree angle, proximal bend with suction control hole	25	74-2083
Oesophageal/ Laryngeal 4mm x 350mm, 30 degree angle, proximal bend with suction control hole	25	74-2087
Oesophageal/ Laryngeal 4mm x 490mm, 30 degree angle, proximal bend with suction control hole	25	74-2090

2mm, 3mm, 4mm Suction Tube Range Specification Sheet**Material Specification**

Product Component	Specification
Tube	Stainless Steel
Moulded handle	Polypropylene (PP)
Pouch Packaging	Tyvek/Film
Outer Box/Carton Packaging	500 Micron White Boxboard

Intended Use

The NETWORK ENT® Suction Sets are designed specifically for ENT procedures and are intended for the removal of excess fluids and debris from the surgical site.

Sterilisation

Products are sterilised by Ethylene Oxide (EO) according to a validated cycle and the requirements of BS EN ISO 11135-1 (current).

Instructions for Use

The instructions for use and suggested surgical technique are supplied in the form of a multi-lingual leaflet with instructions given in diagram form where appropriate. The international symbols used on the packaging are explained in each language. A leaflet is supplied with each product.

Conformity to the European Directives

The Network range of Aspiration Suction Tubes and Fine Ends are identified as surgically invasive devices for transient use and are therefore classified as a Class IIa devices according to Annex IX, Rule 6, of the Medical Devices Directive 93/42/EEC as amended by 2007/47/EC.