



DTR Medical®

AN
INNOVIA
MEDICAL
COMPANY

Sterile Single-use
The surgeon's choice





About Innovia Medical

Innovia Medical is a family of specialty surgical companies that are passionately committed to partnering with medical professionals to help elevate the delivery of patient care and improve clinical outcomes.

ELEVATING CARE. IMPROVING OUTCOMES.

Innovia Medical companies develop and manufacture a portfolio of product solutions designed with both the medical professional and their patients in mind, to deliver care that:

- Saves time for both the medical professional and patient
- Delivers a consistent and predictable care experience
- Improves patient confidence with reduced recovery times and less discomfort

Innovia Medical does this through a commitment to innovation leadership, unsurpassed product quality and a dedication to provide responsive, personalized customer service.

INNOVATION LEADERSHIP.

Innovation leadership is fueled through our dedication to collaboration with medical professionals. We regularly seek feedback and input from physicians, nurses, and other medical professionals through product research, testing and development.

This allows us to truly understand the solutions medical professionals need, and guides our innovative culture focused on inventing new, industry-leading solutions and reinventing current solutions that make ordinary products extraordinary, across the specialty surgery spectrum.

UNSURPASSED PRODUCT QUALITY.

Our products deliver proven solutions and dependability which is the result of rigorous quality control and a relentless commitment to ongoing manufacturing improvements. We use only the highest quality materials and subject all products to rigorous inspection protocols. This is why we're proud to have so many customers who have continued to rely on our solutions for more than forty years.

RESPONSIVE, PERSONALIZED CUSTOMER SERVICE.

From our sales representatives to our customer service team, every associate at Innovia Medical strives to deliver an unrivaled customer experience. This begins by ensuring that each customer receives the product that is right for their specific needs. We provide assistance through technical support and live customer service representatives with extensive product knowledge, available throughout the day. We deliver quick and thorough responses, every time—from order and shipping confirmations to delivery date communications and more ensuring continuous satisfaction.

We know you're committed to providing the best patient care experience possible and that's why we commit to being your trusted partner helping you every step of the way. At Innovia Medical, we are not just products, we are part of the patient care solution.

Learn more at www.innoviamedical.com

Contents

Ear Specula <i>Including Oval, Rosen, Shea & Tumarkin</i>	4-5
Clearway Suction Handle (for foreign body removal)	6
Tongue Depressors <i>Including Lack & Bruenings</i>	7
Nasal Specula <i>Including Hartmann & Thudichum</i>	7-8
Mirrors <i>Including Dental, Laryngeal & Post Nasal</i>	8
Crocodile Forceps & Scissors <i>Including Fine Jaw Micro Crocodile Forceps, Crocodile Micro Forceps & Crocodile Scissors</i>	9-11
Forceps	9-19
Iris Scissors	20
Scissors <i>Including Kilner, Iris, Reynolds, Baby Metzenbaum, Stevens Tenotomy, Wire Cutting, & Mayo Scissors</i>	21-23
Clamps <i>Including Backhaus Towel & Tubing Clamps</i>	24
Tracheal Dilators	25
Hooks & Points <i>Including Cawthorne, Gillies Skin Hook, Rosen Needle & Schuknecht Pick</i>	25-27
Probes <i>Including Jobson Horne & Formby Probes</i>	27
Curettes & Directors <i>Including House & Beckmann Adenoid Curette & Grooved Director</i>	28-29
Retractors <i>Including Maxillofacial, Alms, Kilner "Catspaw" & Baby Roux Retractor</i>	30-31
Rhinology & Otology	32
SNAP, Endoscope Guide	33
Splints	34
Special Suction Handles & Dissectors <i>Including Aspirating, Watson Cheyne & Negus Aspirating Dissector</i>	34-35
Suction Handles <i>Including Barron, Bellucci, Micro Laryngeal, Frazier, Long, Magill, House, Fenestrated, Micro Suction, Suction Regulator, Schuknecht, Sinus, Metal "Glass-Type", Tapered, Yankauer & Zoellner (with Fine Ends) Suction Handles</i>	36-43
Cream Applicators	43
Endoscopic Ear Range <i>Including Curved Zoellner Suction Fine Ends & Curved Endoscopic Suction Tubes</i>	44
Kerrison Rongeurs	45
Scope Warmer	46
Instrument Wipes	46

Ear Specula

- Sterile and non-sterile options available in a range of sizes
- Matt black finish improves light transmission and visibility
- Smooth edges for patient comfort and safety
- Large sizes (wide and short) suitable for giving injections
- Robust Single-use Specula

EAR SPECULA - OVAL 100 (NON-STERILE)

OSP1002	Oval Ear Specula size 2 (3.5mm)	100 units
OSP1003	Oval Ear Specula size 3 (4mm)	100 units
OSP1004	Oval Ear Specula size 4 (6mm)	100 units

EAR SPECULA - OVAL 50

OSP0502	Oval Ear Specula size 2 (3.5mm)	50 units
OSP0503	Oval Ear Specula size 3 (4mm)	50 units
OSP0504	Oval Ear Specula size 4 (6mm)	50 units

EAR SPECULA - OVAL 20

OSP0202	Oval Ear Specula size 2 (3.5mm)	20 units
OSP0203	Oval Ear Specula size 3 (4mm)	20 units
OSP0204	Oval Ear Specula size 4 (6mm)	20 units



EAR SPECULA - ROSEN 100 (NON-STERILE)

ESP1003	Rosen Round Ear Specula 3mm	100 units
ESP1004	Rosen Round Ear Specula 4mm	100 units
ESP1005	Rosen Round Ear Specula 5mm	100 units
ESP1006	Rosen Round Ear Specula 6mm	100 units
ESP1007	Rosen Round Ear Specula 7mm	100 units

EAR SPECULA - ROSEN 50

ESP0503	Rosen Round Ear Specula 3mm	50 units
ESP0504	Rosen Round Ear Specula 4mm	50 units
ESP0505	Rosen Round Ear Specula 5mm	50 units
ESP0506	Rosen Round Ear Specula 6mm	50 units
ESP0507	Rosen Round Ear Specula 7mm	50 units

EAR SPECULA - ROSEN 20

ESP0203	Rosen Round Ear Specula 3mm	20 units
ESP0204	Rosen Round Ear Specula 4mm	20 units
ESP0205	Rosen Round Ear Specula 5mm	20 units
ESP0206	Rosen Round Ear Specula 6mm	20 units
ESP0207	Rosen Round Ear Specula 7mm	20 units



Ear Specula

EAR SPECULA - SHEA 100 (NON-STERILE)

SSP1001	Shea Ear Specula size 1 (3.5 x 4.1mm)	100 units
SSP1002	Shea Ear Specula size 2 (4.5 x 5.5mm)	100 units
SSP1003	Shea Ear Specula size 3 (5.5 x 7.5mm)	100 units
SSP1004	Shea Ear Specula size 4 (6.5 x 8.5mm)	100 units
SSP1005	Shea Ear Specula size 5 (7.5 x 9.5mm)	100 units

EAR SPECULA - SHEA 50

SSP0501	Shea Ear Specula size 1 (3.5 x 4.1mm)	50 units
SSP0502	Shea Ear Specula size 2 (4.5 x 5.5mm)	50 units
SSP0503	Shea Ear Specula size 3 (5.5 x 7.5mm)	50 units
SSP0504	Shea Ear Specula size 4 (6.5 x 8.5mm)	50 units
SSP0505	Shea Ear Specula size 5 (7.5 x 9.5mm)	50 units

EAR SPECULA - SHEA 20

SSP0202	Shea Ear Specula size 2 (4.5 x 5.5mm)	20 units
SSP0203	Shea Ear Specula size 3 (5.5 x 7.5mm)	20 units
SSP0204	Shea Ear Specula size 4 (6.5 x 8.5mm)	20 units
SSP0205	Shea Ear Specula size 5 (7.5 x 9.5mm)	20 units



EAR SPECULA - TUMARKIN 100 (NON-STERILE)

TSP10000	Tumarkin Ear Specula size 0 (3.5mm)	100 units
TSP10005	Tumarkin Ear Specula size 0.5 (4.5mm)	100 units
TSP10010	Tumarkin Ear Specula size 1 (5.5mm)	100 units
TSP10020	Tumarkin Ear Specula size 2 (6mm)	100 units
TSP10030	Tumarkin Ear Specula size 3 (6.5mm)	100 units
TSP10040	Tumarkin Ear Specula size 4 (7mm)	100 units

EAR SPECULA - TUMARKIN 50

TSP05000	Tumarkin Ear Specula size 0 (3.5mm)	50 units
TSP05005	Tumarkin Ear Specula size 0.5 (4.5mm)	50 units
TSP05010	Tumarkin Ear Specula size 1 (5.5mm)	50 units
TSP05020	Tumarkin Ear Specula size 2 (6mm)	50 units
TSP05030	Tumarkin Ear Specula size 3 (6.5mm)	50 units
TSP05040	Tumarkin Ear Specula size 4 (7mm)	50 units

EAR SPECULA - TUMARKIN 20

TSP02000	Tumarkin Ear Specula size 0 (3.5mm)	20 units
TSP02005	Tumarkin Ear Specula size 0.5 (4.5mm)	20 units
TSP02010	Tumarkin Ear Specula size 1 (5.5mm)	20 units
TSP02020	Tumarkin Ear Specula size 2 (6mm)	20 units
TSP02030	Tumarkin Ear Specula size 3 (6.5mm)	20 units
TSP02040	Tumarkin Ear Specula size 4 (7.0mm)	20 units



Clearway Suction Handle

Foreign body removal

- Removes difficult foreign objects from ear and nose
- Effectively latches onto the surface of uneven objects and removes them, such as plugs of dead skin and wax
- Ideal for A&E/ER and MIU
- Can remove the need for ENT surgeon intervention or theatre time
- UK manufactured
- 100% quality inspected
- Our innovative tip, has a soft malleable suctioning cap that forms an air-tight seal to securely grasp foreign objects
- 142° angled handle for improved visibility
- Suction control, particularly useful for smooth, rounded objects
- Universal suction attachment
- Readily available Bubble Tubing or Yankauer Tubing to be cut down at a suitable point once the handle has been removed

“After 60 minutes of trying to retrieve a small ball bearing from a child’s ear, we decided to use the Clearway Suction Handle which removed it in seconds! We now stock these as a matter of course and have used on many occasions since, with excellent feedback from our ENT surgeons”.


East Kent Hospitals University NHS Foundation Trust, 2013



See our video demonstration of the Clearway at:
dtrmedical.com/product/clearway-suction-handle

FOREIGN BODY REMOVAL

SVT1001-S Clearway Suction Handle 10 units



Training poster available on request

Tongue Depressors

LACK TONGUE DEPRESSORS

LAC1001S	Lack Tongue Depressor, Small	20 units
LAC1001M	Lack Tongue Depressor, Medium	20 units
LAC1001L	Lack Tongue Depressor, Large	20 units

- Strong, robust alternative that requires less pressure providing you with an effective substitute over wooden alternatives
- Smooth, matt surface edges allow for misting to assess nasal airflow and minimise patient trauma



BRUENINGS TONGUE DEPRESSOR

BRU2001	Bruenings Tongue Depressor 190mm	20 units
---------	----------------------------------	----------

- Robust alternative to wood
- Matt surface finish prevents misting




Nasal Specula


HARTMANN NASAL SPECULUM

HNS2011-D	Hartmann Nasal Specula, Small	20 units
HNS2021-D	Hartmann Nasal Specula, Medium	20 units
HNS2031-D	Hartmann Nasal Specula, Large	20 units


- Strong stainless steel construction with smooth edges



Small



Medium



Large



Nasal Specula

THUDICHUM NASAL SPECULUM

TNS2003	Thudichum Nasal Specula, multipack (sizes 1, 2 & 3)	10 x 3 units
TNS2000	Thudichum Nasal Specula size 0 (Extra Small)	20 units
TNS2011	Thudichum Nasal Specula size 1 (Small)	20 units
TNS2022	Thudichum Nasal Specula size 2 (Medium)	20 units
TNS2033	Thudichum Nasal Specula size 3 (Large)	20 units

- Balanced sprung blades for ease of use
- Size 0 for paediatric use
- Featuring a strong flexible Stainless Steel construction with smooth edges
- Balanced sprung blades for ease of use
- Multipack containing three individual sizes
- Particularly suitable for Clinic use



Mirrors

FULL RANGE

- Single-use Mirrors with a robust metal construction with smooth edges and secure mirrors

DENTAL MIRRORS

DMM2010	Dental Mirror 10mm diameter	20 units
DMM2015	Dental Mirror 15mm diameter	20 units
DMM2020	Dental Mirror 20mm diameter	20 units
DMM2025	Dental Mirror 25mm diameter	20 units
DMM2030	Dental Mirror 30mm diameter	20 units

LARYNGEAL MIRRORS

LYM2010	Laryngeal Mirror 10mm diameter	20 units
LYM2015	Laryngeal Mirror 15mm diameter	20 units
LYM2020	Laryngeal Mirror 20mm diameter	20 units
LYM2025	Laryngeal Mirror 25mm diameter	20 units
LYM2030	Laryngeal Mirror 30mm diameter	20 units
LYM5015	Laryngeal Mirror 15mm diameter	50 units
LYM5020	Laryngeal Mirror 20mm diameter	50 units
LYM5025	Laryngeal Mirror 25mm diameter	50 units
LYM5030	Laryngeal Mirror 30mm diameter	50 units

POST NASAL MIRRORS

PNM2010	Post Nasal Mirror 10mm diameter	20 units
PNM2012	Post Nasal Mirror 12mm diameter	20 units
PNM2014	Post Nasal Mirror 14mm diameter	20 units



You may also be interested in the ClearView Scope Warmer see p46 for more details

Crocodile Fine Jaw Micro Forceps

Fine specification

- Enables an effective removal of Foreign Bodies from the Ear Cavity
- Commonly used for Wax Removal, Foreign Body Removal and Grommet Placement
- Precise 0.9 x 5mm jaw action that can be used with the smallest Specula
- For Paediatric use or where good vision is critical
- Smooth, matt black finish prevents any reflection when used with a Microscope providing you with enhanced clarity during procedures
- Optimal first-time use and accurate grip
- 100% inspection and dimensional checks on every Instrument

“Really happy with the Fine Jaw
Crocodile Micro Forceps for
Outpatients.”

*ENT Department, NHS Greater Glasgow
and Clyde*



POINT CAMERA TO
VIEW OUR “PRO” FILE

CROCODILE FINE JAW MICRO FORCEPS

CFBF2001-S Crocodile Fine Jaw Micro Forceps black finish

20 units



Crocodile Forceps

CROCODILE MICRO FORCEPS


CF2001-S	Crocodile Micro Forceps steel finish	20 units
CFB2001-S	Crocodile Micro Forceps black finish	20 units

5MM JAW

- Single-use avoids costly repairs to damaged or worn jaws
- Available in steel and matt black finish
- Eliminates risk of contamination from occluded surfaces
- High quality stainless steel with precise jaw action




POINT CAMERA TO
VIEW INFOGRAPHIC



5 x 6mm
jaw



Crocodile Micro Forceps black finish



5 x 6mm
jaw



Crocodile Micro Forceps steel finish


Crocodile Scissors

CROCODILE SCISSORS


CSS1001-S	Crocodile Scissors Straight stainless steel	10 units
CSC1001-S	Crocodile Scissors Curved stainless steel	10 units

STRAIGHT AND CURVED

- Precise cutting action



Curved



Straight




Forceps

HENCKEL AURAL FORCEPS

HAF1001-D	Henckel Aural Forceps	10 units
-----------	-----------------------	----------


- Commonly used for the retrieval and removal of nasal tissue and cartilage
- Fenestrated cup jaws to aid wax removal
- Eliminates risk of contamination from difficult to clean jaw



ORMEROD FORCEPS

OFB1001-S	Ormerod Forceps black finish	10 units
-----------	------------------------------	----------

- High quality stainless steel with precise jaw action
- Matt black finish to minimise reflection
- Eliminates risk of contamination from difficult to clean jaw



Oral Biopsy Forceps

- Ensures a precise biopsy is taken to assist in the diagnosis of Cancer or any other abnormality in the Oral and Laryngeal Mucosa
- With the option 90mm and 120mm, versatility is provided to suit your surgical requirements
- Two cup sizes allowing you flexibility to take a better sample
- Elongated shaft to aid procedure

"We find the Oral Biopsy Forceps very useful for biopsying abnormalities of the tonsil under local anaesthetic in the outpatient clinic, and so avoiding a hospital admission and a visit to the operating theatre."

ENT Department, Doncaster and Bassetlaw Teaching Hospitals NHS Foundation Trust



POINT CAMERA TO
VIEW OUR "PRO" FILE



90mm



ORAL BIOPSY FORCEPS

OBF1001-S	Oral Biopsy Forceps 90mm stainless steel	10 units
OBFL1001-S	Oral Biopsy Forceps Large 120mm stainless steel	10 units

Forceps

SPENCER WELLS ARTERY FORCEPS

SWF2010	Spencer Wells Artery Forceps 120mm Straight	20 units
SWF2015	Spencer Wells Artery Forceps 150mm Straight	20 units
SWF2018	Spencer Wells Artery Forceps 180mm Straight	20 units
SWF2020	Spencer Wells Artery Forceps 200mm Straight	20 units

- Used during surgery to compress the artery, clamp, and seal small blood vessels or hold the artery out of the way
- Gripping jaws and finely serrated blades





Spencer Wells Artery Forceps 200mm Straight

CRILE ARTERY FORCEPS

CFS2001-D	Crile Artery Forceps 140mm Straight	20 units
CFC2001-D	Crile Artery Forceps 140mm Curved on Flat	20 units

- Commonly used for clamping blood vessels or tissue
- Effective locking ratchet and robust grip



Crile Artery Forceps Straight






Crile Artery Forceps Curved

Forceps

MOSQUITO FORCEPS

MQC2012	Micro Mosquito Forceps 120mm Curved	20 units
MQC2025-D	Mosquito Forceps Curved 125mm	20 units
MQC2125-D	Mosquito Forceps Curved Toothed 1x2 125mm	20 units
MQT5001-S	Mosquito Forceps 125mm Straight	50 units

- Frequently used for small blood vessels
- Type of Hemostatic Forcep commonly used to control bleeding
- Robust locking ratchet and serrated jaws for holding Tissue in place



Micro Mosquito Curved 120mm




Mosquito Forceps Straight 125mm

MAGILL CATHETER FORCEPS

MCF11601	Magill Catheter Forceps 160mm	10 units
MCF12001	Magill Catheter Forceps 200mm	10 units
MCF12501	Magill Catheter Forceps 250mm	10 units

- Curved shaft to bend around the tracheal
- Range of lengths available to suit patient needs
- MCF11601 – 2.5mm inner hole diameter x 6mm head width
- MCF12001 – 3.8mm inner hole diameter x 8mm head width
- MCF12501 – 6mm inner hole diameter x 11.5mm head width





Magill Catheter Forceps 250mm

Forceps

KOCHER ARTERY FORCEPS

KAF2014-D Kocher Artery Forceps 140mm 20 units


- Effective locking ratchet and robust grip
- Type of Hemostatic Forceps commonly used to control Bleeding
- 1x2 teeth at the tip to ensure a secure grasp on heavy Tissue and when clamping off bleeding Blood Vessels



KOCHER-OSHSNER FORCEPS

KOF2181 Kocher-Ochsner 1x2 Toothed Forceps Straight 185mm 20 units

- Secure grasp on heavy tissue and when clamping off bleeding Blood Vessels




Keele Mucosectomy Forceps

- Ideal for use during Trans Oral Mucosectomy procedure
- Forceps are used to hold, control, manipulate and stabilise tissue in the Oral Cavity and manoeuvre the tongue during Cauterisation
- Used for removing Carcinomas at the back of the Tongue
- 140° angled tip for enhanced access
- 1x2 Teeth to maintain firm specimen retention

“The Keele Mucosectomy Forceps are my instrument of choice for Endoscopic Mucosectomy. The bespoke designed length and angle allows for excellent manipulation of the tongue and adequate exposure.”



Ajith George, FRCS, Consultant ENT Surgeon, Keele Anatomy and Surgical Training Centre



POINT CAMERA TO
VIEW VIDEO

KEELE MUCOSECTOMY FORCEPS



KMF1025 Keele Mucosectomy Forceps, 250mm, 140° angle 10 units



Forceps

FINE SUTURE FORCEPS

SPF2001-S	Fine Suture Forceps Plain 105mm	20 units
SSF2001-S	Fine Suture Forceps Serrated 105mm	20 units







ADSON FORCEPS

ADS5001-S	Adson Non-Toothed Forceps 120mm	50 units
ATS5001-S	Adson Toothed Forceps 120mm	50 units
ATP2012-D	Adson Forceps 1x2 with plateau 120mm	20 units
ATP2015-D	Adson Forceps 1x2 with plateau 150mm	20 units
ATS2012F	Fine Tip Adson Toothed Forceps 1x2 with plateau 120mm	20 units

- Narrow tips to effectively grip and hold delicate tissue in place
- Flexibility across a wide range of procedures




We also supply Adson Bipolar Forceps, check out our General Surgery catalogue



McINDOE FORCEPS

MDF5001-D	McIndoe Dissecting Forceps 150mm, Non-Toothed	50 units
-----------	---	----------

- Fine, non-toothed serrated jaws are ideal when re-positioning delicate Tissue
- Can be used to hold wound dressings and packing




Forceps

DEBAKEY FORCEPS

DAF2016-D DeBakey Atraumatic Dissecting Forceps 160mm 20 units


- Used to grasp soft tissue and organs during surgery
- Specially designed teeth maintain the grip on the tissues, but do not tear



BABY OVERHOLT FORCEPS

BOC2014 Baby Overholt Artery Forceps Curved 140mm 20 units


- Curved to allow easy access to veins



DENTAL FORCEPS

CDF2001 College Dental Forceps 160mm 20 units


- Precision instrument with firm grip
- Good spring action



GILLIES FORCEPS

GDF2001-D Gillies Toothed Dissecting Forceps 160mm 20 units

- Available with fine, narrow jaws
- 1:2 serrated teeth
- Less pressure required to adequately grip harder material such as skin




Forceps

ALLIS TISSUE FORCEPS

ATF2341 Allis Tissue Forceps 3:4 teeth 160mm 20 units

- Enables a strong grasp in surgical procedures
- Finger ring, ratchet forceps
- Ensure an effective hold of slippery tissue



TILLEY PACKING FORCEPS

TAF2001-S Tilley Aural Forceps 54mm 20 units

TNF2001-S Tilley Nasal Dressing Forceps 76mm 20 units

- Strong slender jaws with good grip and tips that always meet
- Single-use avoids damage or difficult to clean box joints
- Particularly suitable for clinic use




TAF2001-S



WILDE AURAL PACKING FORCEPS

WAF2001-S Wilde Aural Packing Forceps 150mm 20 units

- Used for packing and dressing for ear and aural procedures



EXTENDED RANGE

Iris Scissors

Iris Scissors

- Enhanced accuracy when it comes to delicate surgeries
- First time sharpness which cut even at the tips
- Narrows to a fine and sharp points

SuperCut Iris Scissors

- One serrated and one sharp cutting blade allows a better grip and cut than traditional alternatives
- Sharp tips make the Scissors ideal for precise cutting of soft tissue, sutures, and wound dressings
- Available with straight or curved blades


IRIS SCISSORS

ISS2111-S	Iris Scissors Straight 90mm	20 units
ISC2111-S	Iris Scissors Curved 90mm	20 units
ISS2110-D	Iris Scissors Straight 110mm	20 units
ISC2110-D	Iris Scissors Curved 110mm	20 units


NEW

SUPERCUT IRIS SCISSORS


SISS2111	SuperCut Iris Scissors Straight 90mm	20 units
SISC2111	SuperCut Iris Scissors Curved 90mm	20 units



Iris Scissors Straight 90mm



Curved



Straight



SuperCut Iris Scissors Straight 90mm






IMAGE COMING SOON

Scissors

KILNER SCISSORS

KSS2001-S Kilner Scissors 115mm 20 units



- Stainless steel sharp/sharp scissors
- First time sharpness which cut even at the tips
- Ensure ease when cutting and dissecting delicate tissue in surgical procedures



REYNOLDS SCISSORS

RSC2014-D Reynolds Scissors Curved 145mm 20 units


- Fine tipped scissors enable a precise action when cutting and dissecting soft tissue
- 145mm length helps keep the clinician's hand out of view in demanding procedures



BABY METZENBAUM SCISSORS

ZSC2121-D Baby Metzenbaum Scissors Curved 120mm 20 units

- First time sharpness which cut even at the tips



NEW

SUPERCUT BABY METZENBAUM SCISSORS

SZSC2121 SuperCut Baby Metzenbaum Scissors 120mm, Curved 20 units

- Ideal for fine cutting of soft tissue and blunt dissection
- One serrated blade and one sharp cutting blade allows a superior grip and cut
- Available with straight or curved blades



IMAGE COMING SOON

Scissors

STEVENS TENOTOMY SCISSORS

TSS2001-S	Stevens Tenotomy Scissors 115mm Straight	20 units
TSC2001-S	Stevens Tenotomy Scissors 115mm Curved	20 units

- Smooth and elegant long handle
- Allow the surgeon increased grip and control when cutting in small or narrow areas



Straight



Curved

NEW

SUPERCUT STEVENS TENOTOMY SCISSORS

STSS2001	SuperCut Stevens Tenotomy Scissors 115mm Straight	20 units
STSC2001	SuperCut Stevens Tenotomy Scissors 115mm Curved	20 units

- Ideal for delicate cutting and dissecting soft tissue
- One serrated blade and one sharp cutting blade allows a superior grip and cut
- Available with straight or curved blades





IMAGE COMING SOON

WIRE CUTTING SCISSORS

WCSH2001-S Wire Cutting Scissors Hardened 125mm 20 units


- Improved Scissors enable the cutting of a wide variety of wire
- Used in jaw fixation
- Allows patients or relatives to cut the wire if any vomiting or choking are a risk
- Ideal for cutting standard Cerclage wire efficiently along the whole length of the blade




MAYO SCISSORS

MSS2141	Mayo Scissors Straight Bevelled Blade 140mm	20 units
MSC2141	Mayo Scissors Curved Bevelled Blade 140mm	20 units
MSS2171	Mayo Scissors Straight Bevelled Blade 170mm	20 units
MSC2171	Mayo Scissors Curved Bevelled Blade 170mm	20 units

- Stainless steel sharp/sharp scissors
- First time sharpness which cut even at the tips
- Choice of curved or straight bevelled blades



Straight



Curved





Clamps



CLAMPS

BTC2013-D	Backhaus Towel Clamp 135mm	20 units
TCF2051	Tubing Clamps 150mm	20 units

- Ensures durability and are manufactured using the highest quality stainless steel
- Secure grip for clamping tubing without damage



Backhaus Towel Clamp




Tubing Clamps

Tracheal Dilators


TROUSSEAU TRACHEAL DILATORS

TTD2101	Trousseau Tracheal Dilator 105mm	20 units
TTD2141	Trousseau Tracheal Dilator 140mm	20 units

- Smooth mechanism to ensure high quality performance
- Adult and paediatric sizes



TTD2101



TTD2141



Trousseau Tracheal Dilator

Hooks and Points

CAWTHORNE

CTH2901-D	Cawthorne Hook 1.5mm 90 degree black	20 units
CTP2001-D	Cawthorne Straight Point black	20 units

- Choice of 90 degree hook or point



Hook



Point

Hooks and Points

GILLIES SKIN HOOKS

GHS2001	Gillies Skin Hook Sharp 182mm	20 units
GHB2001	Gillies Skin Hook Semi-Sharp 182mm	20 units

- Available in both sharp and Semi-sharp offering flexibility in a variety of procedures
- The sharp hook enables enhanced accuracy, grip and precision
- The semi-sharp hook is suitable for holding skin whilst taking a biopsy, minimising patient trauma



ROSEN NEEDLES

RSN2016	Rosen Needle 160mm	20 units
RSNF2016	Rosen Needle Fine 160mm	20 units

- Fine variant to complement existing Rosen Needle for more precise procedures
- Constructed to create a clean incision in the Ear Drum when used in Surgeries such as Tympanoplasty
- Ideal for the removal of Wax
- Sharp tip you can rely on in delicate Operations
- Developed and produced using Consultant feedback




Hooks and Points

SCHUKNECHT PICK

SPS2001 Schuknecht Pick Straight 160mm 20 units

- Round knurled handle for grip
- Shaft angled 18 degrees



Probes

JOBSON HORNE PROBES

JHPS5001-S	Jobson Horne Probe 127mm (5 inch)	50 units
JHPL5001-S	Jobson Horne Probe 178mm (7inch)	50 units
JHPS10001-S	Jobson Horne Probe 127mm (5 inch)	100 units
JHPL10001-S	Jobson Horne Probe 178mm (7inch)	100 units

- Easy to bend to shape to match anatomy and improve vision
- Fine working ends for use in even the smallest ear specula
- Jobson double ended loop and cotton wool grip
- Single-use avoids need to re-straighten or risk breakage during use



Jobson Horne Probe 127mm



Jobson Horne Probe 178mm

FORMBY CERUMEN HOOK & PROBE

FCH5001-S Formby Cerumen Hook and Probe 178mm 50 units

- Used to remove foreign objects from the ear
- Most commonly used for the removal of earwax



House Curette

- Slender arms for easy access and good visibility within the ear cavity
- Scoop tips enable a quick and efficient removal of tough cortical bone
- Sharp double ended 2/3mm scoops gently scrape at the bone during middle ear procedures, ensuring a clean cut
- Flat handles with grip allow for superior control and effective scraping
- Used in Stapedotomy to improve hearing in patients
- Suitable for Endoscopic Ear Surgery

“The curettes are very good. I use them in endoscopic ear surgery (EES). The main difference apart from disposable & always sharp is the fact that length & size is suitable for EES. Because of the flat middle bit, we can generate more torque when curetting.”

*Mr Arunachalam Iyer,
Monklands Hospital*



HOUSE CURETTE

HCD22030 House Curette, Double Ended 2.0/3.0mm Scoops, 180mm

20 units

Curettes

BECKMANN ADENOID CURETTE

BAC2010	Beckmann Adenoid Curette, 10mm Curette Edge, 220mm	20 units
BAC2013	Beckmann Adenoid Curette, 13mm Curette Edge, 220mm	20 units
BAC2015	Beckmann Adenoid Curette, 15mm Curette Edge, 220mm	20 units
BAC2017	Beckmann Adenoid Curette, 17mm Curette Edge, 220mm	20 units

- Sharp and precise scoops to enable an effective adenoidectomy
- Serrated handle to provide a secure grip whilst scraping tissue




Directors

GROOVED DIRECTOR

GDR2001	Grooved Director 130mm	20 units
---------	------------------------	----------

- Provides a clear line of sight





Retractors

MAXILLOFACIAL RETRACTOR

MFR0410 Maxillofacial Retractors x 4 with silicone loop 10 units

- Stainless steel dual hook with blunt points
- Latex-free 2mm silicone 180mm ties
- Pack of 4 for wound size flexibility



ALMS RETRACTOR

ALR1101 Alms Retractor Large 103mm 10 units



- Rigid stainless steel construction
- High quality four hook proximal screw action that locks the arms into place



KILNER "CATSPAW" RETRACTOR

KCP2001 Kilner "Catspaw" Retractor 155mm 20 units

- First time sharpness for precision
- Precise 2.2mm width prong end and hook blade guarantees reliability and performance



Retractors

BABY ROUX RETRACTOR

BRR2013 Baby Roux Retractor 130mm 20 units

- Curved, atraumatic ends enable you to effectively retract tissue with minimal trauma, as well as allowing greater accessibility in surgical procedures




Rhinology and Otology

MYRINGOTOME NARROW ARROW & MICROEDGE BLADE HANDLES

MNA2001	Myringotome Narrow Arrow	20 units
BL-1425	MicroEdge Blade Handle, reusable, non-sterile	1 unit

- Ability to bend to required angle
- Individually packed to avoid cross contamination
- First-time sharpness for precision



Myringotome Narrow Arrow



Myringotome Handle

TILLEY LICHTWITZ SINUS PUNCH

TLP1001 Tilley Lichtwitz Sinus Punch 10 units


- First time sharpness for precision



VENTILATION TUBES AND GROMMETS

VT-0201-01	Shepard vent tubes with wire 1.14mm lumen fluoroplastic	6 units
VT-0202-01	Shepard vent tubes 1.14mm lumen fluoroplastic	6 units
VT-0203-01	Shepard vent tubes with wire 1.0mm lumen fluoroplastic	6 units
VT-0204-01	Shepard vent tubes 1.0mm lumen fluoroplastic	6 units
VT-0802-01	Shah with tab 1.14mm lumen fluoroplastic	6 units
VT-0804-01	Shah 1.14mm lumen fluoroplastic	6 units
VT-0400-01	T-tube 12 x 10 x 1.14mm lumen silicone	6 units
VT-0401-01	T-tube 5 x 8 x 1.14mm lumen silicone	6 units
VT-0403-01	T-tube Grommet 1.27mm lumen silicone	6 units
VT-0300-01	Paparella 1 1.14mm lumen silicone	6 units
VT-0501-01	Armstrong bevelled straight shank 7.5 x 1.14mm lumen fluoroplastic	6 units
VT-0503-01	Armstrong Grommet, 1.14mm lumen fluoroplastic	6 units
VT-0503-50	Armstrong Grommet, 1.14mm lumen fluoroplastic	50 units

- Comprehensive range of designs, sizes and materials
- Available with or without wires



SNAP, Endoscope Guide

- Protects from airborne risks and diseases
- A one-way valve that can be retrofitted to any fluid resistant Surgical Mask
- Designed with the aim to keep Clinicians safe and so better able to deliver high quality care to their patients
- Facilitates the movement of a Nasendoscope up to 4mm in diameter
- Ensures safe passage of a Nasendoscopy through a Surgical Mask, creating an aperture for the Endoscope to pass through into the Nasal Cavity
- Reduces the risks associated with the patient sneezing
- Scientifically proven to reduce the prevalence of AGP (Aerosol Generating Particles) during examination
- Provides better safety and care for both the Patient & Clinician
- Clinical Otolaryngology Peer Review Paper and Trial
- Filtrex Optical Particle Report for SNAP
- MET 1952-2591 Splash Resistance Test Report
- SNAP White Paper August 2020
- UK Manufactured
- Winner of the 2020 MediWales "Partnership with the NHS" Award
- Winner of the 2021 "Special Focus On COVID" Meridian Award
- Continuing ENT Diagnosis Throughout COVID-19 Safely webinar
- UK Patent Application Number 2105626.2
- Exclusively distributed by DTR Medical



"On reflection we should have been masking patients for Endoscopy many years ago"

Ajith George, FRCS, Consultant ENT Surgeon

SNAP, ENDOSCOPE GUIDE

SNAP40 SNAP, Endoscope Guide 40 Unit



POINT CAMERA TO
VIEW VIDEO

Splints

NASAL SPLINTS

10-1500-05KL Denver Nasal Splint Kit Large

5 units

- Non-stick splint reduces adhesion
- Slit improves insertion and removal
- Pre-punched holes provided to ease suturing



Special Suction Handles and Dissectors

FREER ASPIRATING DISSECTOR

FDS0501 Freer Aspirating "Hockey Stick" Dissector

5 units

- Freer Aspirating "Hockey Stick" for nasal procedures
- Optimum blade sharpness for precision work
- Removes risk from very difficult to clean suction blades



[View Our Negus Aspirating Dissector p35](#)

WATSON CHEYNE DISSECTOR


WCD2225 Watson Cheney Dissector, 225mm

20 units

- Flexible device can be used for the mobilisation and removal of degenerative tissue
- Very fine blunt ends with slender curves which allow for high precision and control
- Avoids the use of traditional re-usable alternatives that are "knobbly and fat"



Pointed end
0.7mm width



Flat end
1.7mm width



EXTENDED RANGE

Negus Aspirating Dissectors

- Optimum blade sharpness for precision work
- Removes risk from very difficult to clean suction blades
- Enlarged hole size increases the rear aperture size
- Front edge of hole is 0.75mm closer to the blade tip
- Enables an effective dissection and manipulation of Nasal Mucosa with suction of fluids



Good suction flow at the blade tip giving better visibility



Less prone to blockage and subsequent theatre delays



Controlled suction for delicate tissue



NDS0511
2.25mm hole
diameter



NDS0521
2.25mm hole
diameter



Negus One-Hole

Negus Two-Hole

NEGUS ASPIRATING DISSECTORS

NDS0511	Negus Aspirating Dissector with One Hole	5 units
NDS0521	Negus Aspirating Needle with Two Holes	5 units

Suction Handles

BARRON SUCTION HANDLES

BSH2003	Barron Suction Handle 1mm (3Fg)	20 units
BSH2004	Barron Suction Handle 1.35mm (4Fg)	20 units
BSH2005	Barron Suction Handle 1.67mm (5Fg)	20 units
BSH2006	Barron Suction Handle 2mm (6Fg)	20 units
BSH2007	Barron Suction Handle 2.3mm (7Fg)	20 units

- Complete range of sizes from 3 to 7Fg with suction control
- Working length of 77mm



MICRO LARYNGEAL SUCTION HANDLES

240mm		
BLSF2240	Micro Laryngeal Suction Handle 240 x 2mm OD (6Fg)	20 units
280mm		
BLS2280	Micro Laryngeal Suction Handle 280 x 2mm OD (6Fg)	20 units
330mm		
BLS2330	Micro Laryngeal Suction Handle 330 x 2mm OD (6Fg)	20 units
BLSF2339	Micro Laryngeal Suction Handle 330 x 3mm OD (9Fg)	20 units
360mm		
BLSF2369	Micro Laryngeal Suction Handle, 360 x 3mm OD (9Fg)	20 units

- Effective alternative to reprocessing very long narrow lumen
- Range of sizes from 240mm (6Fg) to 360mm (9Fg)




EXTENDED RANGE

Frazier Suction Handles

- Strong construction suitable for bending to shape during procedure
- Single-use avoids the need to restraighten or risk breakage
- Tip options include non-fenestrated, fenestrated with two holes and cross fenestrated with four holes
- Expanded range with 150° degree bend
- Pre-bent handles
- 5mm (15Fg) variant now available

Supplied straight* and pre bent
Bend to shape for your procedure

**Excluding 12Fg and 15Fg range*



FRAZIER SUCTION HANDLES

2mm (6Fg)

FSH2006	Frazier Suction Handle 200 x 2mm OD (6Fg) straight	20 units
FSH2006-B	Frazier Suction Handle 190 x 2mm OD (6Fg) 150° bend	20 units
FSF2006	Frazier Suction Handle 200 x 2mm OD (6Fg) straight, fenestrated	20 units
FSF2046	Frazier Suction Handle 200 x 2mm OD (6Fg) straight, cross fenestrated	20 units
FSN2426	Frazier Suction Handle 184 x 2mm OD (6Fg), 42° bend, without suction control	20 units

2.66mm (8Fg)

FSH2008	Frazier Suction Handle 200 x 2.66mm OD (8Fg) straight	20 units
FSH2008-B	Frazier Suction Handle 190 x 2.66mm OD (8Fg) 150° bend	20 units
FSF2008	Frazier Suction Handle 200 x 2.66mm OD (8Fg) straight, fenestrated	20 units
FSF2048	Frazier Suction Handle 200 x 2.66mm OD (8Fg) straight, cross fenestrated	20 units

3mm (9Fg)


FSH2009	Frazier Suction Handle 200 x 3mm OD (9Fg) straight	20 units
FSH2009-B	Frazier Suction Handle 190 x 3mm OD (9Fg) 150° bend	20 units
FSF2009	Frazier Suction Handle 200 x 3mm OD (9Fg) straight, fenestrated	20 units
FSF2049	Frazier Suction Handle 200 x 3mm OD (9Fg) straight, cross fenestrated	20 units
FSN2429	Frazier Suction Handle 184 x 3mm OD (9Fg), 42° bend, without suction control	20 units

4mm (12Fg)

FSH2012	Frazier Suction Handle 200 x 4mm OD (12Fg), 150° bend	20 units
FSF2012	Frazier Suction Handle 200 x 4mm OD (12Fg), fenestrated, 150° bend	20 units
FSF2412	Frazier Suction Handle 200 x 4mm OD (12Fg), cross fenestrated, 150° bend	20 units
FSN24212	Frazier Suction Handle 184 x 4mm OD (12Fg), 42° bend, without suction control	20 units

5mm (15Fg)

FSH2015	Frazier Suction Handle 200 X 5mm OD (15Fg), 150° bend	20 units
---------	---	----------



Suction Handles



LONG SUCTION HANDLES

450mm		
BLS2450	Bronchoscopy Suction Handle 450 x 4mm OD (12Fg)	20 units
550mm		
BLS2550	Oesophageal Suction Handle 550 x 4mm OD (12Fg)	20 units
560mm		
SGH2001	Sigmoidoscope Suction Handle, 560 x 6mm OD (18Fg)	20 units


- Effective alternative to reprocessing very long narrow lumen



Bronchoscopy Suction Handle 450 x 4mm



Sigmoidoscope Suction Handle, 560 x 6mm OD (18Fg)




Oesophageal Suction Handle 550 x 4mm

MAGILL SUCTION HANDLES

MGH2010	Magill Suction Handle 1mm Lumen 210 x 3mm OD	20 units
MGH2015	Magill Suction Handle 1.5mm Lumen 210 x 3.5mm OD	20 units
MGH2020	Magill Suction Handle 2mm Lumen 210 x 4mm OD	20 units

- Atraumatic tip



Magill Suction Handle 1mm Lumen 210 x 3mm OD

Suction Handles

HOUSE SUCTION HANDLES

HST0118/18	House Suction and Irrigation Handle (1.8 x 1.8mm)	20 units
HST0125/20	House Suction and Irrigation Handle (2.5 x 2mm)	20 units
HST0130/20	House Suction and Irrigation Handle (3 x 2mm)	20 units
HST0145/20	House Suction and Irrigation Handle (4.5 x 2mm)	20 units
HST0145/30	House Suction and Irrigation Handle (4.5 x 3mm)	20 units
HST0150/45	House Suction and Irrigation Handle (5 x 4.5mm)	20 units


- Dual purpose suction and irrigation with suction control
- Complete range of sizes from 1.8 to 5mm
- Single-use overcomes repairs to tips damaged through drilling
- Particularly suitable for collection of bone dust



FENESTRATED HOUSE SUCTION HANDLE

HSFL118/18	House Suction Handle Fenestrated (1.8 x 1.8mm)	20 units
------------	--	----------

- 1.8 x 1.8mm fenestration for further suction



Suction Handles

EXTENDED RANGE

MICRO SUCTION TUBES

86mm

MSV508113	Micro Suction (Verhoeven) Red 86 x 2.3mm OD (13g)	50 units
MSV508113R	Micro Suction (Verhoeven) Red 86 x 2.3mm OD (13g) with Suction Regulator	50 units
MSV508114	Micro Suction (Verhoeven) Green 86 x 2.1mm OD (14g)	50 units
MSV508114R	Micro Suction (Verhoeven) Green 86 x 2.1mm OD (14g) with Suction Regulator	50 units
MSV508116	Micro Suction (Verhoeven) Cream 86 x 1.65mm OD (16g)	50 units
MSV508117	Micro Suction (Verhoeven) White 86 x 1.47mm OD (17g)	50 units
MSV508119	Micro Suction (Verhoeven) Brown 86 x 1.05mm OD (19g)	50 units

110mm

MSV50114	Micro Suction (Verhoeven) Green 110 x 2.1mm OD (14g)	50 units
MSV50117	Micro Suction (Verhoeven) White 110 x 1.47mm OD (17g)	50 units
MSV50119	Micro Suction (Verhoeven) Brown 110 x 1.05mm OD (19g)	50 units
MSV50122	Micro Suction (Verhoeven) Black 110 x 0.7mm OD (22g)	50 units

120mm

MSV50214	Micro Suction (Verhoeven) Green 120 x 2.1mm OD (14g)	50 units
----------	--	----------

- Stainless steel tube with colour coded luer connector for microsuction use
- Also suitable for use with syringe as irrigation device
- Available in 86, 110 and 120mm lengths



UNIVERSAL ADAPTER

SUCTION REGULATOR

DSRS5001	Suction Regulator, Stainless Steel with Improved Grip	50 units
----------	---	----------

- House style regulator with control
- Luer connector suitable for Verhoeven and other micro suction tubes



For Sterile Tube Sets, see our General Surgery catalogue




Suction Handles

SCHUKNECHT SUCTION HANDLES

SHK2020	Schuknecht Suction Handle 20swg	20 units
SHK2024	Schuknecht Suction Handle 24swg	20 units

- Single-use removes need to clean a very narrow lumen
- Working length of 77mm




FULL RANGE

SINUS SUCTION HANDLES

SSH2606	Sinus Suction Handle (6Fg) 65 degrees	20 units
SSH2609	Sinus Suction Handle (9Fg) 65 degrees	20 units
SSH2909	Sinus Suction Handle (9Fg) 90 degrees	20 units
SSH2612	Sinus Suction Handle (12Fg) 65 degrees	20 units
SSH2912	Sinus Suction Handle (12Fg) 90 degrees	20 units



65 degrees



5.5mm
tip length


90 degrees



METAL "GLASS-TYPE" SUCTION HANDLE

GSH2001	Metal "Glass-Type" Suction handle 200 x 5mm OD (15Fg)	20 units
---------	---	----------

- Metal version of traditional glass suction handle
- Suitable for where wide bore is required



22mm tip
length



Suction Handles

TAPERED SUCTION HANDLES

190mm

TDH2195	Tapered Suction Handle 190mm size 5 with tear drop control	20 units
TDH2196	Tapered Suction Handle 190mm size 6 with tear drop control	20 units
TDH2197	Tapered Suction Handle 190mm size 7 with tear drop control	20 units

220mm

TDH2225	Tapered Suction Handle 220mm size 5 with tear drop control	20 units
TDH2226	Tapered Suction Handle 220mm size 6 with tear drop control	20 units
TDH2227	Tapered Suction Handle 220mm size 7 with tear drop control	20 units

- Tear drop control for adjustable suction




CHOICE OF ROSES

YANKAUER SUCTION HANDLES

YST1000	Yankauer Suction Handle 305mm with rose	10 units
YST1000R	Yankauer Suction Handle 305mm with modified rose	10 units

- Ideal when plastic handles are inappropriate
- Detachable rose with three (YST1000) or six holed modified rose (YST1000R)
- Ideal for tonsillectomy




Suction Handles

ZOELLNER SUCTION HANDLES

SMS4400-S	Zoellner Suction Handle	50 units
SMS4418-S	Zoellner Suction Handle with 18swg Fine End	50 units

- Available with or without fine ends
- Ergonomic balanced design for ease of use
- Tube rotation suits left or right handed users




Zoellner Suction Handle without Fine End

ZOELLNER FINE ENDS

SM4018	Zoellner Fine Ends size 18 swg	50 units
SM4020	Zoellner Fine Ends size 20 swg	50 units
SM4022	Zoellner Fine Ends size 22 swg	50 units
SM4024	Zoellner Fine Ends size 24 swg	50 units
SM4026	Zoellner Fine Ends size 26 swg	50 units

- Robust metal taper suitable for refitting during procedure
- Extra long taper gives a good fit with any suction handle
- Range of sizes from 18 to 26 standard wire gauge



Cream Applicators

CREAM APPLICATORS


ECA5001-S	Cream Applicator 86mm long x 6fg	50 units
-----------	----------------------------------	----------

- Cream or ointment applicator for use with luer lock syringe




Endoscopic Ear Range

- Extensive range for growth in ear endoscopy which provide clear visibility at the procedure site
- Overcomes cleaning difficulty and cross contamination with narrow lumen curved tubes



Curved Zoellner Suction Fine End

Curved Zoellner Suction Fine End



Curved Endoscopic Suction Handles

CURVED ZOELLNER SUCTION FINE ENDS

SMES4483-D	Curved Endoscopic Zoellner Fine Ends, 3mm Bend Radius, 18 swg with Zoellner Suction Handle	20 units
SMES4486-D	Curved Endoscopic Zoellner Fine Ends, 6mm Bend Radius, 18 swg with Zoellner Suction Handle	20 units

CURVED ENDOSCOPIC SUCTION TUBES

SME4183	Curved Endoscopic Zoellner Suction Fine End, 3mm Bend Radius, 18 swg	20 units
---------	--	----------

Kerrison Rongeur


- Latest 100mm variant with 90-degree up-bite complements our existing range
- Well balanced stainless steel construction
- Specifically designed for cutting and removing small pieces of bone during surgery
- Eliminates the risk of cross contamination between procedures



KRUF11029



KRUF11039



KERRISON RONGEURS

KRU1101	Kerrison Rongeur Punch 1mm wide, length 200mm with 40 degree up bite	10 units
KRU1201	Kerrison Rongeur Punch 2mm wide, length 200mm with 40 degree up bite	10 units
KRU1301	Kerrison Rongeur Punch 3mm wide, length 200mm with 40 degree up bite	10 units
KRU1401	Kerrison Rongeur Punch 4mm wide, length 200mm with 40 degree up bite	10 units
KRU1501	Kerrison Rongeur Punch 5mm wide, length 200mm with 40 degree up bite	10 units
KRUF11029	Kerrison Rongeur Punch 2mm wide, length 100mm with 90 degree up bite	10 units
KRUF11039	Kerrison Rongeur Punch 3mm wide, length 100mm with 90 degree up bite	10 units

Scope Warmer

THEATRE EFFICIENCY

SCOPE WARMER

31SW10	ClearView Anti-Fog Scope Warmer, sterile packed for Single-use	10 units
31SW30	ClearView Anti-Fog Scope Warmer, sterile packed for Single-use	30 units

- Simple preheating to eliminate condensation
- Secures the surgeon's vision
- Warms full scope length to 37°C within 90 seconds, staying above body temperature for 120 minutes
- Long-acting to avoid disruption or repeated defogging




Instructional poster available
Contact us on +44 (0)1792 797910

Instrument Wipes

NETCELL® PVA INSTRUMENT WIPE

35-900	PVA Instrument Wipe	20 units
--------	---------------------	----------

- Quickly removes protein, blood and particles without damage to the instrument
- Manufactured from ultra-smooth, biocompatible PVA sponge that will not snag
- Soft, strong, flexible, fibre and snag free material



"Supporting visibility, quality and safety in your supply chain"***

Our instruments' packaging features the latest machine readable GS1-standard barcoded labels for product traceability.

The newly incorporated labelling technology is a timely development with an increasing number of health providers worldwide adopting an 'eProcurement' strategy and will provide machine readable product information on all labels to meet expected EU and other regulators demands.

Associated benefits include:

- Improve patient care and reduce costs
- Improve stock control
- One global language
- Accurately track equipment

"Unique device identification (UDI) compliance represents a significant development for our customers. Automation of information flow will have influences across the whole purchasing and supply process, providing the benefits of superior instrument traceability and improved patient safety."

"We are working closely with Kodit UDI software to provide our labelling solutions to meet the latest international regulations and we serve as a reference site for other manufacturers to learn from."

DTR Medical

For more information on our new product labelling, please contact Customer Services on +44 (0) 1792 797910 or dtrinfo@innoviamedical.com

** - GS1 – Global Traceability Standards

Matching reusable standards

We benchmark against reusable British and international standards:

BS 5194-2:1989

Specification for Instruments with Pivot Joints

BS 5194-3:1985

Specification for Dissecting Forceps

BS 5194-4:1985

Specification for Scissors, Shears and other Jointed Instruments

We operate our own UK production facilities equipped with ISO class 7 cleanrooms

DTR Medical

ENT - Single-Use Instruments

DTR Medical Ltd

An Innovia Medical Company

17 Clarion Court, Enterprise Park, Swansea SA6 8RF UK

Telephone: +44 (0)1792797910

dtrinfo@innoviamedical.com

US 1.888.229.2875 | **International** +44 (0)1792797910

www.innoviamedical.com/single-use-instruments