




Product Overview

The NETWORK ENT® Fine Ends range has been specifically designed to work alongside the NETWORK ENT® range of suction tubes to provide increased precision during delicate procedures. These bulk non-sterile products are designed for inclusion in custom procedure packs.

- Provides pin point suction for micro ENT procedures
- Manufactured with soft coloured-coded seal, for easy identification and superior fit
- Compatible with other makes of suction tube
- Supplied non-sterile, single use only
- Manufactured in the UK

Fine End Size Options

Product	Size	Unit of Sale	Code
Fine End (Purple)	18G / 1.0 mm	Bag of 100	74-2020NS
Fine End (Orange)	20G / 0.8 mm	Bag of 100	74-2021NS
Fine End (Blue)	22G / 0.65 mm	Bag of 100	74-2022NS
Fine End (Green)	24G / 0.5 mm	Bag of 100	74-2023NS

Fine Ends Specification Sheet
BULK KITPACK PRODUCT**Material Specification****Product Component****Specification**

Tube	Stainless Steel
Moulded handle	Polypropylene (PP)

Intended Use

Network Aspirating and Suction tubes are designed particularly for ENT applications and when attached to a vacuum suction system, are used to remove excess fluids and debris from the surgical site. The Vari-Tube Fine ends are inserted into the tube ends of the Vari-Tube Suction handle to reduce the diameter of the tube for precision and pinpoint suctioning during delicate otological procedures.

Sterilisation

Products are supplied non-sterile, and are suitable for sterilisation by EO gas, gamma irradiation and e-beam. Devices are designated for Single Use Only.

Instructions for Use

No instructions for use are provided.

Conformity to the European Directives

The Network range of Aspiration Suction Tubes and Fine Ends are identified as surgically invasive devices for transient use and are therefore classified as a Class IIa devices according to Annex IX, Rule 6, of the Medical Devices Directive 93/42/EEC as amended by 2007/47/EC.