GENERAL & VASCULAR SURGERY



Sterile Single-use Enhancing surgical workflow





About Innovia Medical

Innovia Medical is a family of specialty surgical companies that are passionately committed to partnering with medical professionals to help elevate the delivery of patient care and improve clinical outcomes.

ELEVATING CARE. IMPROVING OUTCOMES.

Innovia Medical companies develop and manufacture a portfolio of product solutions designed with both the medical professional and their patients in mind, to deliver care that:

- Saves time for both the medical professional and patient
- Delivers a consistent and predictable care experience
- Improves patient confidence with reduced recovery times and less discomfort

Innovia Medical does this through a commitment to innovation leadership, unsurpassed product quality and a dedication to provide responsive, personalized customer service.

INNOVATION LEADERSHIP.

Innovation leadership is fueled through our dedication to collaboration with medical professionals. We regularly seek feedback and input from physicians, nurses, and other medical professionals through product research, testing and development.

This allows us to truly understand the solutions medical professionals need, and guides our innovative culture focused on inventing new, industry-leading solutions and reinventing current solutions that make ordinary products extraordinary, across the specialty surgery spectrum.

UNSURPASSED PRODUCT QUALITY.

Our products deliver proven solutions and dependability which is the result of rigorous quality control and a relentless commitment to ongoing manufacturing improvements. We use only the highest quality materials and subject all products to rigorous inspection protocols. This is why we're proud to have so many customers who have continued to rely on our solutions for more than forty years.

RESPONSIVE, PERSONALIZED CUSTOMER SERVICE.

From our sales representatives to our customer service team, every associate at Innovia Medical strives to deliver an unrivaled customer experience. This begins by ensuring that each customer receives the product that is right for their specific needs. We provide assistance through technical support and live customer service representatives with extensive product knowledge, available throughout the day. We deliver quick and thorough responses, every time—from order and shipping confirmations to delivery date communications and more ensuring continuous satisfaction.

We know you're committed to providing the best patient care experience possible and that's why we commit to being your trusted partner helping you every step of the way. At Innovia Medical, we are not just products, we are part of the patient care solution.

Learn more at www.innoviamedical.com

Contents

Forceps Crile Artery, Mosquito, Kocher Artery, Kocher-Ochsner, Adson, Mcindoe, Fine Suture, Debakey, Gillies, Baby Overholt, Allis Tissue, Rampley Sponge Holders, Spencer Wells & Strassburg Forceps	4-7
Needle Holders Mayo-Hegar, Crilewood, Halsey, Ryder, Kilner & DeBakey Needle Holders	8-9
Scissors Kilner, Iris, Vannas, Mayo, Reynolds, Metzenbaum, Stevens Tenotomy & SuperCut Scissors	9-14
Clamps Backhaus Towel & Tubing Clamps	15
Bipolar Forceps Standard, Adson, Cawthorne, Bipolar Cable, Monopolar Forceps & Bayonet Bipolar Forceps	16-18
Curettes Beckmann Adenoid Curette	19
Clamp Covers Sterile & Non-Sterile Clamp Covers	19
Silicone Slings	20
General & Vascular Heparin Flushing Needles	20
Tibbs Arterial Cannula	21
Phlebectomy Hooks	22
Dissectors Watson Cheyne Dissector	22
Retractors Kilner Alae Hook Retractor	22
Suction Handles Bellucci, Long, Micro Laryngeal, Broncoscopy, Oesophageal, Sigmoidoscope, Frazier, Micro, Suction Regulator, Poole & Yankauer	23-20
Tip Shields	26
Tubing Products Silicone & Sterile Silicone Tube Sets	27-30
Others Silicone Bands, Blue Clamps & Scalpel Blade Handles	31
Scope Warmer	32
Instrument Wipes	32

CRILE ARTERY FORCEPS CFS2001-D Crile Artery Forceps 145mm Straight 20 units CFC2001-D Crile Artery Forceps 145mm Curved on Flat 20 units

- Commonly used for clamping blood vessels or tissue
- Effective locking ratchet and robust grip








MOSQUITO FORCEPS			
MQT5001-S	Mosquito Forceps 125mm Straight	50 units	
MQC2012	Micro Mosquito Forceps 120mm Curved	20 units	
MQC2025-D	Mosquito Forceps Curved 125mm	20 units	
M0C2125-D	Mosquito Forcens Curved Toothed 1x2 125mm	20 units	

- Frequently used for small blood vessels
- Type of Hemostatic Forceps commonly used to control bleeding
- Robust locking ratchet and serrated jaws for holding tissue in place








KOCHER ARTERY FORCEPS

KAF2014-D Kocher Artery Forceps 140mm 20 units


- Effective locking ratchet and robust grip
- Type of Hemostatic Forceps commonly used to control Bleeding
- 1×2 teeth at the tip, our Forceps ensure a secure grasp on heavy Tissue and when clamping off bleeding Blood Vessels





KOCHER-OCHSNER FORCEPS

KOF2181 Kocher-Ochsner 1x2 Toothed Forceps Straight 185mm 20 units







ADSON FO	RCEPS	
ADS5001-S	Adson Non-Toothed Forceps 120mm	50 units
ATS5001-S	Adson Toothed Forceps 120mm	50 units
ATP2012-D	Adson Forceps 1x2 with plateau 120mm	20 units
ATP2015-D	Adson Forceps 1x2 with plateau 150mm	20 units
ATS2012F	Fine Tip Adson Toothed Forceps 1x2 with plateau 120mm	20 units

- Narrow tips to effectively grip and hold delicate tissue in place
- Flexibility across a wide range of procedures





Adson Toothed Forceps



Adson Non-Toothed Forceps



Adson Forceps with Plateau



Fine Tip Adson Toothed Forceps

McINDOE FORCEPS

MDF5001-D McIndoe Dissecting Forceps 150mm, Non-Toothed 50 units


- Fine, non-toothed serrated jaws are ideal when re-positioning delicate Tissue
- Can be used to hold wound dressings and packing





FINE SUTURE FORCEPS SPF2001-S Fine Suture Forceps Plain 105mm 20 units SSF2001-S Fine Suture Forceps Serrated 105mm 20 units

- Provides flexibility required in delicate procedures
- Deliver an excellent grip to help secure tissue and suture materials







Fine Suture Forceps Serrated

DEBAKEY FORCEPS

DAF2016-D DeBakey Atraumatic Dissecting Forceps 160mm 20 units

- Used to grasp soft tissue and organs during surgery
- Specially designed teeth maintain the grip on the tissues, but do not tear







GILLIES FORCEPS

GDF2001-D Gillies Toothed Dissecting Forceps 160mm 20 units

- Available with fine, narrow jaws
- 1:2 serrated teeth
- Less pressure required to adequately grip harder material such as skin








BABY OVERHOLT FORCEPS

BOC2014 Baby Overholt Artery Forceps Curved, 140mm 20 units

■ Curved to allow easy access to veins






ALLIS TISSUE FORCEPS

ATF2341 Allis Tissue Forceps 3:4 teeth, 160mm 20 units

- Enables a strong grasp in surgical procedures
- Finger ring, ratchet forceps
- Ensure an effective hold of slippery tissue





RAMPLEY SPONGE HOLDERS

RSH2181 Rampley Sponge Holder 180mm 20 units RSH2241 Rampley Sponge Holder 240mm 20 units

 Strong jaws with an effective locking ratchet for effectively holding swabs and gripping tissue during procedures

■ Robust grip





SPENCER WELLS FORCEPS

SWF2010	Spencer Wells Artery Forceps 120mm Straight	20 units
SWF2015	Spencer Wells Artery Forceps 150mm Straight	20 units
SWF2018	Spencer Wells Artery Forceps 180mm Straight	20 units
SWF2020	Spencer Wells Artery Forceps 200mm Straight	20 units

- Used during surgery to compress the Artery, clamp, and seal small Blood Vessels or hold the Artery out of the way
- Gripping jaws and finely serrated blades





Needle Holders

MAYO-HE	GAR NEEDLE HOLDERS	
MNH2013	Mayo-Hegar Needle Holder 130mm	20 units
MNH2015	Mayo-Hegar Needle Holder 150mm	20 units
MNHG2013-D	Mayo-Hegar Needle Holder with groove 130mm	20 units

- Engineered with serrated jaws that allow for an effective grip
- Central groove to ensure the Needle can be securely held




MNH2015




MNHG2013-D

NEEDLE HO	OLDERS	
CNH2015	Crilewood Needle Holder 150mm	20 units
HNH2013	Halsey Needle Holder 130mm	20 units
RNH2015-D	Ryder Needle Holder 150mm	20 units

- Ensure a firm grip when Suturing
- Robust locking ratchet







Needle Holders

NEEDLE HOLDERS KNH2013 Kilner Needle Holder 135mm 20 units DNH2018-D DeBakey Needle Holder 180mm 20 units

Designed with narrow, atraumatic teeth to ensure a secure grasp on fine Sutures and delicate tissue whilst minimising damage





Kilner Needle Holder 135mm



Scissors

KILNER SCISSORS KSS2001-S Kilner Scissors 115mm 20 units

- Stainless steel sharp/sharp scissors
- First time sharpness which cut even at the tips
- Ensure ease when cutting and dissecting delicate tissue in surgical procedures





Scissors

IRIS SCISS	SORS	
ISS2111-S	Iris Scissors Straight 90mm	20 units
ISC2111-S	Iris Scissors Curved 90mm	20 units
ISS2110-D	Iris Scissors Straight 110mm	20 units
ISC2110-D	Iris Scissors Curved 110mm	20 units

- Enhanced accuracy when it comes to delicate surgeries
- First time sharpness which cut even at the tips
- Narrows to a fine and sharp point







NEW

SUPERCU	JT IRIS SCISSORS	
SISS2111	SuperCut Iris Scissors Straight 90mm	20 units
SISC2111	SuperCut Iris Scissors Curved 90mm	20 units

- One serrated and one sharp cutting blade allows a better grip and cut than traditional alternatives
- Sharp tips make the Scissors ideal for precise cutting of soft tissue, sutures, and wound dressings
- Available with straight or curved blades





IMAGE COMING SOON


Vannas Scissors


- High quality stainless steel cross action scissors with fine sharp tips
- Precise blade action
- Curved option offers greater flexibility for delicate procedures
- Ideal for your most delicate practice, such as the removal of fine Sutures
- Equipped with a ridged handle and a double leaf spring mechanism for self-opening action, ensuring enhanced control

Straight









Precise blade

VANNAS SCISSORS


VSS2081	Vannas Scissors Straight 85mm	20 units
VSS2121	Vannas Scissors Straight 125mm	20 units
VSC2081	Vannas Scissors Curved 85mm	20 units

MAYO SC	ISSORS	
MSS2141	Mayo Scissors Straight Bevelled Blade 140mm	20 units
MSC2141	Mayo Scissors Curved Bevelled Blade 140mm	20 units
MSS2171	Mayo Scissors Straight Bevelled Blade 170mm	20 units
MSC2171	Mayo Scissors Curved Bevelled Blade 170mm	20 units

- Stainless steel sharp/sharp scissors
- First time sharpness which cut even at the tips
- Choice of curved or straight bevelled blades



Straight







REYNOLDS SCISSORS

RSC2014-D Reynolds Scissors Curved 145mm 20 units

- Fine tipped scissors enable a precise action when cutting and dissecting soft tissue
- 145mm length helps keep the clinician's hand out of view in demanding procedures








Metzenbaum Scissors

- First time performance for cutting delicate tissue and blunt dissection
- Specially engineered to be fine and precise



NEW

SUPERCUT BABY METZENBAUM SCISSORS SZCS2121 SuperCut Baby Metzenbaum Scissors 120mm, Curved 20 units

Ideal for fine cutting of soft tissue and blunt dissection

■ One serrated blade and one sharp cutting blade allows a superior grip and cut

IMAGE COMING SOON

STEVENS TENOTOMY SCISSORS


TSS2001-S Stevens Tenotomy Scissors 115mm Straight 20 units TSC2001-S Stevens Tenotomy Scissors 115mm Curved 20 units

- Stainless steel sharp/sharp Scissors
- First time sharpness which cut even at the tips
- Smooth and elegant long handle
- Allows increased grip and control when cutting in small or narrow areas











NEW

SUPERCUT STEVENS TENOTOMY SCISSORS

STSS2001 SuperCut Stevens Tenotomy Scissors 115mm Straight 20 units STSC2001 SuperCut Stevens Tenotomy Scissors 115mm Curved 20 units

- Ideal for delicate cutting and dissecting soft tissue
- One serrated blade and one sharp cutting blade allows a superior grip and cut
- Available with straight or curved blades






IMAGE COMING SOON


Clamps


CLAMPS		
BTC2013-D	Backhaus Towel Clamp 135mm	20 units
TCF2051	Tubing Clamps 150mm	20 units

- Ensures durability and are manufactured using the highest quality stainless steel
- Secure grip for clamping tubing without damage










Bipolar Forceps -

STANDARD E	BIPOLAR FORCEPS		
SU770513044	Standard Bipolar Forceps 130 x 0.5mm, 3m, 4mm Plugs	20 units	
SU770515044	Standard Bipolar Forceps 150 x 0.5mm, 3m, 4mm Plugs	20 units	
SU770520044	Standard Bipolar Forceps 200 x 0.5mm, 3m, 4mm Plugs	20 units	
SU771015044	Standard Bipolar Forceps 150 x 1mm, 3m, 4mm Plugs	20 units	
SU771520044	Standard Bipolar Forceps 200 x 1.5mm, 3m, 4mm Plugs	20 units	
SU772015044	Standard Bipolar Forceps 150 x 2mm, 3m, 4mm Plugs	20 units	
SU772020044	Standard Bipolar Forceps 200 x 2mm, 3m, 4mm Plugs	20 units	

- Fine stainless steel tips with a range of sizes from 0.5mm 2mm
- Adjustable and universal 4mm plug that is tested to comply with the latest 'Flying Leads' UK legislation
- High quality cables and connectors deliver consistently repeatable results





ADSON BIPO	OLAR FORCEPS	
SU790515044	Adson Bipolar Forceps 150 x 0.5mm, 3m, 4mm Plugs	20 units
SI 1791015044	Adson Rinolar Forcens 150 x 1mm 3m 4mm Plugs	20 units





Bipolar Forceps

CAWTHORNE BIPOLAR FORCEPS SU760508044 Cawthorne Bipolar Forceps 80 x 0.5mm, 3m, 4mm Plugs 20 units SU761008044 Cawthorne Bipolar Forceps 80 x 1mm, 3m, 4mm Plugs 20 units SU760511044 Cawthorne Bipolar Forceps 110 x 0.5mm, 3m, 4mm Plugs 20 units SU761011044 Cawthorne Bipolar Forceps 110 x 1mm, 3m, 4mm Plugs 20 units



BIPOLAR CABLE

BPC5030 Bipolar Cable, 3.0m lead, 4mm pins UK/EU Fitting 50 units





MONOPOLAR FORCEPS


SU671220044 Monopolar Forceps 200mm Straight x 1.2mm tip serrated 20 units





Bayonet Bipolar Forceps

- Fine stainless steel tips with 1mm or 2mm options
- Adjustable and universal 4mm plug that is tested to comply with the latest 'Flying Leads' UK legislation
- High quality cables and connectors deliver consistently repeatable results



Curettes -

BECKMANN ADENOID CURETTE			
BAC2010	Beckmann Adenoid Curette, 10mm Curette Edge, 220mm	20 units	
BAC2013	Beckmann Adenoid Curette, 13mm Curette Edge, 220mm	20 units	
BAC2015	Beckmann Adenoid Curette, 15mm Curette Edge, 220mm	20 units	
BAC2017	Beckmann Adenoid Curette, 17mm Curette Edge, 220mm	20 units	

- Sharp and precise scoops to enable an effective Adenoidectomy
- Serrated handle to provide a secure grip whilst scraping Tissue





Clamp Covers

STERILE CLAMP COVERS			
CCR2031	Clamp Cover Red (Small) 50mm	20 units	
CCW2041	Clamp Cover White (Med) 120mm	20 units	
CCB2051	Clamp Cover Blue (Large) 120mm	20 units	

CLAMP COVERS (NON-STERILE) SM3040 White - Medium 5 x 1m

Blue - Large

SM3050

Soft silicone jaw covers ensure good grip with reduced

5 x 1m

 Custom sterile packs cut to preferred lengths available on request



Silicone Slings —

SILICONE	SILICONE SLINGS			
1.5mm Slings	1.5mm Slings			
SM2100	Multicoloured Slings 1.5mm packed x 5	20 units		
SM2103	Red Slings 1.5mm	25 units		
SM2104	Blue Slings 1.5mm	25 units		
SM2105	Yellow Slings 1.5mm	25 units		
SM2106	White Slings 1.5mm	25 units		
2mm Slings				
SM2002 SM2008 SM2030	Multicoloured Slings 2mm packed x 5 Sling Multipack 2mm, (Red x 3, Blue x 3, Yellow x 1, White x 1) Red Slings 2mm	20 units 20 units 25 units		
SM2040	Blue Slings 2mm	25 units		
SM2050	Yellow Slings 2mm	25 units		
SM2060	White Slings 2mm	25 units		



■ Radiopaque and biocompatible medical grade silicone

Colour coded for your specific procedure:

Blue: VeinRed: ArteryYellow: NerveWhite: Tendon



General and Vascular ————

	HEPARIN	FLUSHING NEEDLES	
	HFN2001	Heparin Flushing Needle 48mm	20 units
l	HFN2081	Heparin Flushing Needle 85mm	20 units


- Available in two distinct lengths to enable better access to the Vein
- Single-use removes the risk of cross-contamination making it very suitable for your high-risk cases





UNIQUE RANGE **Tibbs Arterial Cannula**

- Unique range of Single-Use Tibbs Arterial Cannula with a high-quality stainless-steel construction
- 8 variants are provided, available in Micro, Small, Large, X Large, and XX Large
- The Cone and Olive variants increase flexibility, ensuring you get the correct size for the specific vessel when flushing them with heparin solution
- UK manufactured



TIBBS STANDARD SIZES

TAC20LC	Tibbs Arterial Cannula Large Cone, 48mm	20 units
TAC20SC	Tibbs Arterial Cannula Small Cone, 48mm	20 units
TAC20L0	Tibbs Arterial Cannula Large Olive, 48mm	20 units
TAC20S0	Tibbs Arterial Cannula Small Olive, 48mm	20 units







TIBBS SPECIAL SIZES

TAC20M0	Tibbs Arterial Cannula Micro Olive 2.5mm OD, 48mm	20 units
TAC20XLC	Tibbs Arterial Cannula X Large Cone 8mm OD extra-long, 69mm	20 units
TAC20XXLC	Tibbs Arterial Cannula XX Large Cone 12.25mm OD extra-long, 69mm	20 units
TAC2003	Tibbs Arterial Cannula Set of 3 sizes (Large Olive, Large Cone & X Large Cone)	20 units



Phlebectomy Hooks

PHLEBEC [*]		
PBH2001	Phlebectomy Hook Size 1 (2mm)	20 units
PBH2002	Phlebectomy Hook Size 2 (1.5mm)	20 units
PBH2003	Phlebectomy Hook Size 3 (1mm)	20 units

- Allows for an effective removal of Varicose Veins that lie underneath the surface of the leg
- Handle shaped for comfort and control during procedures
- Readily available in 1.0, 1.5 and 2.0mm Hook sizes to ensure a precise incision







Dissectors

WATSON CHEYNE DISSECTOR

WCD2225 Watson Cheyne Dissector, 225mm 20 units

- Flexible device can be used for the mobilisation and removal of degenerative Tissue in General Surgery
- Very fine blunt ends with slender curves which allow for high precision and control
- Avoids the use of traditional re-usable alternatives that are "knobbly and fat"



Retractors

KILNER ALAE HOOK RETRACTOR

KHR2001 Kilner Alae Hook Retractor 20 units

■ Effective retraction for surgical exposureon





Suction Handles ————

LONG SUC	LONG SUCTION HANDLES		
240mm			
BLSF2240	Micro Laryngeal Suction Handle 240 x 2mm OD (6Fg)	20 units	
280mm			
BLS2280	Micro Laryngeal Suction Handle 280 x 2mm OD (6Fg)	20 units	
330mm			
BLS2330	Micro Laryngeal Suction Handle 330 x 2mm OD (6Fg)	20 units	
BLSF2339	Micro Laryngeal Suction Handle 330 x 3mm OD (9Fg)	20 units	
360mm			
BLSF2369	Micro Laryngeal Suction Handle, 360 x 3mm 0D (9Fg)	20 units	
450mm			
BLS2450	Broncoscopy Suction Handle 450 x 4mm OD (12Fg)	20 units	
550mm			
BLS2550	Oesophageal Suction Handle 550 x 4mm OD (12Fg)	20 units	
560mm SGH2001	Sigmoidoscope Suction Handle, 560 x 6mm OD (18Fg)	20 units	

- Effective alternative to reprocessing very long narrow lumen
- Range of sizes from 240mm (6Fg) to 560mm (18Fg)



EXTENDED RANGE

Frazier Suction Handles

- Strong construction suitable for bending to shape during procedure
- Single-use avoids the need to restraighten or risk breakage
- Tip options include non-fenestrated, fenestrated with two holes and cross fenestrated with four holes
- Expanded range with 150° degree bend
- Pre-bent handles
- 5mm (15Fg) variant now available

Supplied straight* and pre bent Bend to shape for your procedure *Excluding 12Fg range

With suction control Cross fenestrated With suction control Without suction 150° bend 42° bend

FRAZIER SUCTION HANDLES

2mm (6Fg) FSH2006 FSH2006-B FSF2006 FSF2046 FSN2426	Frazier Suction Handle 200 x 2mm OD (6Fg) straight Frazier Suction Handle 190 x 2mm OD (6Fg) 150° bend Frazier Suction Handle 200 x 2mm OD (6Fg) straight, fenestrated Frazier Suction Handle 200 x 2mm OD (6Fg) straight, cross fenestrated Frazier Suction Handle 184 x 2mm OD (6Fg), 42° bend, without suction control	20 units 20 units 20 units 20 units 20 units
2.66mm (8Fg FSH2008 FSH2008-B FSF2008 FSF2048	Frazier Suction Handle 200 x 2.66mm OD (8Fg) straight Frazier Suction Handle 190 x 2.66mm OD (8Fg) 150° bend Frazier Suction Handle 200 x 2.66mm OD (8Fg) straight, fenestrated Frazier Suction Handle 200 x 2.66mm OD (8Fg) straight, cross fenestrated	20 units 20 units 20 units 20 units
3mm (9Fg) FSH2009 FSH2009-B FSF2009 FSF2049 FSN2429	Frazier Suction Handle 200 x 3mm OD (9Fg) straight Frazier Suction Handle 190 x 3mm OD (9Fg) 150° bend Frazier Suction Handle 200 x 3mm OD (9Fg) straight, fenestrated Frazier Suction Handle 200 x 3mm OD (9Fg) straight, cross fenestrated Frazier Suction Handle 184 x 3mm OD (9Fg), 42° bend, without suction control	20 units 20 units 20 units 20 units 20 units
4mm (12Fg) FSH2012 FSF2012 FSF2412 FSN24212	Frazier Suction Handle 200 x 4mm OD (12Fg), 150° bend Frazier Suction Handle 200 x 4mm OD (12Fg), fenestrated, 150° bend Frazier Suction Handle 200 x 4mm OD (12Fg), cross fenestrated, 150° bend Frazier Suction Handle 184 x 4mm OD (12Fg), 42° bend, without suction control	20 units 20 units 20 units 20 units
5mm (15Fg) FSH2015	Frazier Suction Handle 200 X 5mm OD (15Fg), 150° bend	20 units

Without suction

Suction Handles ————

EXTENDED RANGE

MICRO SUC	MICRO SUCTION TUBES		
86mm			
MSV508113	Micro Suction (Verhoeven) Red 86 x 2.3mm OD (13g)	50 units	
MSV508113R	Micro Suction (Verhoeven) Red 86 x 2.3mm OD (13g) with Suction Regulator	50 units	
MSV508114	Micro Suction (Verhoeven) Green 86 x 2.1mm OD (14g)	50 units	
MSV508114R	Micro Suction (Verhoeven) Green 86 x 2.1mm OD (14g) with Suction Regulator	50 units	
MSV508116	Micro Suction (Verhoeven) Cream 86 x 1.65mm OD (16g)	50 units	
MSV508117	Micro Suction (Verhoeven) White 86 x 1.47mm OD (17g)	50 units	
MSV508119	Micro Suction (Verhoeven) Brown 86 x 1.05mm OD (19g)	50 units	
110mm		,	
MSV50114	Micro Suction (Verhoeven) Green 110 x 2.1mm OD (14g)	50 units	
MSV50117	Micro Suction (Verhoeven) White 110 x 1.47mm OD (17g)	50 units	
MSV50119	Micro Suction (Verhoeven) Brown 110 x 1.05mm OD (19g)	50 units	
MSV50122	Micro Suction (Verhoeven) Black 110 x 0.7mm OD (22g)	50 units	
120mm MSV50214	Micro Suction (Verhoeven) Green 120 x 2.1mm OD (14g)	50 units	



- Also suitable for use with syringe as irrigation device
- Available in 86, 110 and 120mm lengths



UNIVERSAL ADAPTER

SUCTION REGULATOR

DSRS5001 Suction Regulator, Stainless Steel with Improved Grip 50 units

- House style regulator with control
- Luer connector suitable for Verhoeven and other micro suction tubes



Suction Handles -

SPECIAL SIZES

	POOLE SU		
	PSHN1001	Poole Suction Handle Neonatal	10 units
	PSH1001	Poole Suction Handle	10 units
l	PSHL1001	Poole Suction Handle Long	10 units

- Stainless steel for balance and precision
- Supplied with detachable outer to spread suction
- Available in three lengths 190, 240, 330mm to suit patients from neonates to high BMI adults





CHOICE OF ROSES

YANKAUER SUCTION HANDLES						
YST1000	Yankauer Suction Handle 305mm with rose	10 units				
YST1000R	Yankauer Suction Handle 305mm with modified rose	10 units				

- Detachable rose with three (YST1000) or six holed modified rose (YST1000R)
- Ideal for tonsillectomy





Tip Shields

TIP SHIEL	.DS	
SM5020	Flat 13mm width Brown	100 units
SM5050	Round 6mm diameter Yellow	100 units
SM5060	Round 4mm diameter Red	100 units
SM5080	Round 2mm diameter Blue	100 units
SM5090	Round 1mm diameter Blue	100 units

- Supplied non-sterile
- Wide choice of colour, size, and shape
- Cost effective instrument protection
- Protects staff during reprocessing and sterilisation



Tubing Products ————

EASY TO USE SIZE GUIDE

SILICONE	TUBING		
SM135	ID 0.51	OD 0.94	15m roll
SM155	ID 0.64	OD 1.19	15m roll
SM205	ID 1.02	OD 2.16	15m roll
SM265	ID 1.57	OD 2.41	15m roll
SM285	ID 1.57	OD 3.18	15m roll
SM305	ID 1.98	OD 3.18	15m roll
SM335	ID 3.35	OD 4.65	15m roll
SM365	ID 3.18	OD 6.35	15m roll
SM405	ID 4.8	OD 7.9	15m roll
SM445	ID 6.4	OD 9.5	15m roll
SM465	ID 6.4	OD 11.1	15m roll
SM485	ID 6.4	OD 12.7	15m roll
SM525	ID 9.5	OD 12.7	15m roll
SM545	ID 9.5	OD 14.3	15m roll

- Supplied non-sterile
- Supplied in 15 metre packs for cutting to size
- Suitable for custom procedure trays



SUCTION H	SUCTION HANDLE TUBING RECOMMENDATIONS									
Tubing (non-sterile)	Tubing Diameters	Spigot	Sterile DSTS Code	Suitable for use Suction Handles						
SM405	4.8mm ID x 7.9mm OD	4-6mm	DSTS40506	Aspirating Dissectors Bronchoscopy Frazier Magill Microlaryngeal Oesophageal Poole Sinus Suction Regulator						
SM445	6.4mm ID x 9.5mm OD	7-11mm	DSTS44506	All of the above (exc Sinus) plus Barron, Schuknecht and Tapered						

Sterile Single-Use Silicone Tube Sets

HOW TO USE: First 6 digits + order code followed by inner diameter (ID) and outer diameter (OD) then length (mm) e.g. DST205 + 1004 = DST2051004 (1.02 ID x 2.16 OD) 40mm

DSTS = Spigot	Inner Diameter	Outer Diameter																			
CODE	mm	mm	40	75	150	200	250	300	400	450	600	700	750- 760	800	1000	1100	1350- 1450	1500- 1750	2000	2500	3600
DST105	0.30	0.64						03													
DST205	1.02	2.16	10												10						
DST265									54												
DST285	1.57	3.18					1or2		03	04		07									
DST305	1.98	3.18				02															
DST325	2.65	4.88						03													
DST365	3.20	6.40												08		11			20		
DSTS365	3.20	6.40				02															
DST405	4.80	7.90			015			03			06		76-S		10		13	15 or 17			
DST415	5.00	8.30									06							15			
DSTS425	4.80	9.50											76								
DST445	6.40	9.50						03				07			10			15		25	36
DSTS445	6.40	9.50									06				10						
DST465	6.40	11.10								04											36
DST545	9.50	12.70		007																	
DSTS545	9.50	14.30					02														
DST625	12.70	17.50																		25	
PACK SIZE KEY	10 (units)	20 (units)																			

YOUR SINGLE-USE SOLUTION WITH CONFIDENCE

Tubing Products

EXTENDED RANGE

STERILE SI	LICONE TUBE SETS	
DST105 DST10503	Sterile Silicone Tube 0.3 ID x 0.64 OD x 300mm long	10 units
DST205 DST2051004 DST20510	Sterile Silicone DCR Tube 1.02 ID x 2.16 OD x 40mm long Sterile Silicone Tube 1.02 ID x 2.16 OD x 1000mm long	10 units 10 units
DST265 DST26554	Sterile Silicone Tube 1.57 ID x 2.41 OD x 400mm long	10 units
DST285 DST2851 DST2852 DST28503 DST28504 DST28507	Sterile Silicone Tube 1.57 ID x 3.18 OD x 250mm long Sterile Silicone Tube 1.57 ID x 3.18 OD x 250mm long 2 pieces Sterile Silicone Tube 1.57 ID x 3.18 OD x 300mm long Sterile Silicone Tube 1.57 ID x 3.18 OD x 450mm long Sterile Silicone Tube 1.57 ID x 3.18 OD x 700mm long	10 units 10 units 10 units 10 units 10 units
DST305 DST30502	Sterile Silicone Tube 1.98 ID x 3.18 OD x 200mm long	10 units
DST325 DST32503	Sterile Silicone Tube 2.65 ID x 4.88 OD x 300mm long	10 units
DST365 DSTS36502 DST36508 DST36511 DST36520	Sterile Silicone Tube 3.2 ID x 6.4 OD x 200mm long with Spigot Sterile Silicone Tube 3.2 ID x 6.4 OD x 800mm long Sterile Silicone Tube 3.2 ID x 6.4 OD x 1100mm long - House Sterile Silicone Tube 3.2 ID x 6.4 OD x 2000mm long	20 units 10 units 20 units 10 units

- Sterile for extra confidence where lengths over 270mm are impossible to validate
- Cost effective alternative to in-house production
- Full range of tube sizes cut to consistent lengths
- Tailored to suit needs including mixed packs
- Selected tube sets available with Spigot





STERILE TUBING AND SPIGOT STARTER PACKS

Code	Description	Pack size
DSTS40506	Sterile Silicone Tube 4.8mm ID x 7.9mm OD x 600mm long with Spigot	10 units
DSTS44506	Sterile Silicone Tube 6.4mm ID x 9.5mm OD x 600mm long with Spigot	10 units

Tubing Products —

CTEDII E CI	LICONE TUBE SETS (CONTINUED)	
	LICONE TOBE SETS (CONTINUED)	
DST405 DST540506 DST405015 DST40503 DST40506 DST40576-S DST40510 DST40513	Sterile Silicone Tube 4.8 ID x 7.9 OD x 1000mm long Sterile Silicone Tube 4.8 ID x 7.9 OD x 1350mm long	10 units 10 units 10 units 10 units 10 units 20 units
DST40515 DST40517	Sterile Silicone Tube 4.8 ID x 7.9 OD x 1500mm long Sterile Silicone Tube 4.8 ID x 7.9 OD x 1750mm long	20 units 20 units
DST415 DST41506 DST41515	Sterile Silicone Tube 5.0 ID x 8.3 OD x 600mm long Sterile Silicone Tube 5.0 ID x 8.3 OD x 1500mm long	10 units 20 units
DST425 DSTS42576	Sterile Silicone Tube 4.8 ID x 9.5 OD x 760mm long with Spigot	10 units
DST445 DST544506 DST44503 DST44507 DST44510 DST544510 DST44515 DST44525 DST44536	Sterile Silicone Tube 6.4 ID x 9.5 OD x 600mm long with Spigot Sterile Silicone Tube 6.4 ID x 9.5 OD x 300mm long Sterile Silicone Tube 6.4 ID x 9.5 OD x 700mm long Sterile Silicone Tube 6.4 ID x 9.5 OD x 1000mm long Sterile Silicone Tube 6.4 ID x 9.5 OD x 1000mm long with Spigot Sterile Silicone Tube 6.4 ID x 9.5 OD x 1500mm long Sterile Silicone Tube 6.4 ID x 9.5 OD x 2500mm long Sterile Silicone Tube 6.4 ID x 9.5 OD x 2500mm long Sterile Silicone Tube 6.4 ID x 9.5 OD x 3600mm long	10 units 10 units 10 units 10 units 10 units 20 units 20 units
DST465 DST46504 DST46536 DST545 DST545007 DSTS54502	Sterile Silicone Tube 6.4 ID x 11.1 OD x 450mm long Sterile Silicone Tube 6.4 ID x 11.1 OD x 3600mm long Sterile Silicone Tube 9.5 ID x 14.3 OD x 75mm long Sterile Silicone Tube 9.5 ID x 14.3 OD x 250mm long with Spigot	20 units 10 units 10 units 10 units
DST625 DST62525	Sterile Silicone Tube 12.7 ID x 17.5 OD x 2500mm long	20 units

Others

SILICONE BANDS

DSB3211-S Silicone Bands 2.16 x 110mm diameter x 3 20 units

- Sterile for extra confidence
- Cost effective alternative to in-house production
- Latex free to replace elastic bands
- Customised options available



BLUE CLAMPS

DBC2051-SBlue Clamp Mediplast20 x 1 unitDBC5051-SBlue Clamp Mediplast10 x 5 units



UNIVERSAL ADD ON

SCALPEL BLADE HANDLES

SBH2031 Scalpel Blade Handle No 3 20 units
SBH2041 Scalpel Blade Handle No 4 20 units

- Secure grip
- Eliminates cross infection



Scope Warmer

THEATRE EFFICIENCY

SCOPE WARMER

31SW10	ClearView Anti-Fog Scope Warmer, sterile packed for Single-use	10 units
31SW30	ClearView Anti-Fog Scope Warmer, sterile packed for Single-use	30 units

- Simple preheating to eliminate condensation
- Secures the surgeon's vision
- Warms full scope length to 37°C within 90 seconds, staying above body temperature for 120 minutes
- Long-acting to avoid disruption or repeated defogging







Instructional poster available Contact us on +44 (0)1792 797910

Instrument Wipes

NETCELL® PVA INSTRUMENT WIPE

35-900 PVA Instrument Wipe 20 units

- Quickly removes protein, blood and particles without damage to the instrument
- Manufactured from ultra-smooth, biocompatible PVA sponge that will not snag
- Soft, strong, flexible, fibre and snag free material



Labelling

"Supporting visibility, quality and safety in your supply chain"**

Our instruments' packaging features the latest machine readable GS1-standard barcoded labels for product traceability.

The newly incorporated labelling technology is a timely development with an increasing number of health providers worldwide adopting an 'eProcurement' strategy and will provide machine readable product information on all labels to meet expected EU and other regulators demands.

Associated benefits include:

- Improve patient care and reduce costs
- Improve stock control
- One global language
- Accurately track equipment

"Unique device identification (UDI) compliance represents a significant development for our customers. Automation of information flow will have influences across the whole purchasing and supply process, providing the benefits of superior instrument traceability and improved patient safety."

"We are working closely with Kodit UDI software to provide our labelling solutions to meet the latest international regulations and we serve as a reference site for other manufacturers to learn from."

DTR Medical

For more information on our new product labelling, please contact Customer Services on +44 (0) 1792 797910 or dtrinfo@innoviamedical.com

** - GS1 – Global Traceability Standards

Matching reusable standards

We benchmark against reusable British and international standards:

BS 5194-2:1989

Specification for Instruments with Pivot Joints

BS 5194-3:1985

Specification for Dissecting Forceps

BS 5194-4:1985

Specification for Scissors, Shears and other Jointed Instruments

We operate our own UK production facilities equipped with ISO class 7 cleanrooms

DTR Medical

General & Vascular - Single-Use Instruments

DTR Medical Ltd An Innovia Medical Company

17 Clarion Court, Enterprise Park, Swansea SA6 8RF UK Telephone: +44 (0)1792797910

dtrinfo@innoviamedical.com

US 1.888.229.2875 | **International** +44 (0)1792797910 www.innoviamedical.com/single-use-instruments