




DTR Medical®

AN
INNOVIA
MEDICAL
COMPANY

THE 'PRO'FILE

'Needles' to Say, First Time Sharpness is Guaranteed - Try Our Sterile Single-use Aspirating Needle



The Design Checklist

- ✓ Stainless steel cannula with 18g needle for laparoscopic aspiration
- ✓ Extra long 335mm length suitable for all patients of any BMI
- ✓ Available in 3.5 and 5mm diameter cannulas complete with Luer connection

Clinicians' Experiences

Don't just take our word for it ...

"Sharpness is essential"

James Cook University Hospital

"Accurate during procedure, effective with good results"

Mid Cheshire NHS Trust

"Sharp is what we need every time"

Royal Alexandra Paisley Hospital

Clinical Value - Optimum First-time Sharpness of Needle for Ease of Use



TIME

A delicate Single-use instrument that is sharp every time

With prompt delivery guaranteed, we endeavour to meet your service needs. Avoiding costly delays, reducing hospital stock holding and improving lead time for instruments



LIFE

Single-use avoids difficult to clean cannula and needle with stepped joint

Single-use eliminates the risk of cross contamination



COST

This is a cost-effective alternative in comparison to poor quality

No burden of sterilisation costs

No repair costs for bent or damaged needles

No shortage of instruments whilst undergoing repair

YOUR SINGLE-USE SOLUTION WITH CONFIDENCE

t +44 (0)1792 797910 e marketing@dtrmedical.com w www.dtrmedical.com

