



DTR Medical®

AN
INNOVIA
MEDICAL
COMPANY



THE 'PRO'FILE

Ensure a Quick & Efficient Procedure with our Sterile Single-use Audiology Range of Instruments




Ear Specula

Rosen, Oval, Shea & Tumarkin available



**Crocodile Fine Jaw
Micro Forceps**



**Clearway
Suction Handle**

The Design Checklist

- ✓ Plethora of effective Single-use Audiology Instruments that assist in the restoration of hearing through a Hearing Aid when the mechanism becomes sclerosed or fused together
- ✓ Enable an efficient removal of hardened wax, foreign objects, and dead skin cells that attach to the Tympanic Membrane

Clearway Suction Handle

- ✓ Innovative tip has a soft malleable suctioning cap that forms an air-tight seal to securely grasp
- ✓ 142° angled handle for improved visibility
- ✓ Suction control, particularly useful for smooth, rounded objects
- ✓ 360° rotation for improved clinician visibility

Crocodile Fine Jaw Micro Forceps

- ✓ Precise 0.9 x 5mm jaw action that can be used with the smallest Specula
- ✓ For Paediatric use or where good vision is critical
- ✓ Smooth, matt black finish prevents any reflection when used with a Microscope providing you with enhanced clarity during procedures
- ✓ Optimal first-time use and accurate grip

Formby Cerumen Hook & Probe

- ✓ 360° tip for effective grip of cotton wool used in a wide range of Otological procedures
- ✓ 3mm outer diameter loop for cerumen extraction
- ✓ Atraumatic fine tips to minimise patient discomfort
- ✓ Double ended for ease of use

Ear Specula

- ✓ Matt black finish improves light transmission and vision
- ✓ Smooth edges for patient comfort and safety
- ✓ Large sizes (wide and short) suitable for giving injections

Zoellner Suction Handle

- ✓ Available with or without Fine Ends
- ✓ Ergonomic balanced design for ease of use
- ✓ Tube rotation suits left or right-handed users

YOUR SINGLE-USE SOLUTION WITH CONFIDENCE

t +44 (0)1792 797910 **e** dtrmarketing@innoviamedical.com **w** www.dtrmedical.com





DTR Medical®

AN
INNOVIA
MEDICAL
COMPANY



THE 'PRO'FILE

Clinical Value – Allow an Effective Restoration of Hearing



TIME

Avoids unsuccessful attempts of wax, dead skin or foreign body retrieval

Put your time to better use and try our range of Sterile Single-use Audiology Instruments that work first time, every time

Single-use prevents the need for you to attempt to re-straighten the Probes, this both damages them and wastes clinical time



LIFE

More comfortable experience for patient as repeated attempts not only cause further swelling and bleeding but also compromise patient cooperation

Eliminates the risk of cross-contamination



COST

No cost of replacing re-usable Instruments that can be easily damaged

Diminished re-processing, repairs, and replacement

Clinicians' Experiences

Don't just take our word for it ...

"DTR Medical's Ear Specula offer the best size and quality."

Hull University Teaching Hospitals NHS Trust

"After 60 minutes of trying to retrieve a small ball bearing from a child's ear, we decided to use the Clearway Suction Handle which removed it in seconds! We now stock these as a matter of course and have used on many occasions since, with excellent feedback from our ENT surgeons."

East Kent Hospitals University NHS Foundation Trust

YOUR SINGLE-USE SOLUTION WITH CONFIDENCE

t +44 (0)1792 797910 **e** dtrmarketing@innoviamedical.com **w** www.dtrmedical.com

