



DTR Medical®

AN
INNOVIA
MEDICAL
COMPANY

THE **'PRO'**FILE

'Fine'-ally Forceps Suitable for Delicate Ear Procedures - Use our Sterile Single-use Crocodile Fine Jaw Micro Forceps



The Design Checklist

- ✓ The best Single-use option available, making them the product of choice when removing Foreign Bodies from the Ear Cavity
- ✓ Commonly used for Wax Removal, Foreign Body Removal and Grommet Placement
- ✓ Precise 0.9 x 5mm jaw action that can be used with the smallest Specula
- ✓ For Paediatric use or where good vision is critical
- ✓ Smooth, matt black finish prevents any reflection when used with a Microscope providing you with enhanced clarity during procedures
- ✓ Optimal first-time use and accurate grip
- ✓ 100% inspection and dimensional checks on every Instrument

YOUR SINGLE-USE SOLUTION WITH CONFIDENCE

t +44 (0)1792 797910 e marketing@dtrmedical.com w www.dtrmedical.com





DTR Medical®

AN
INNOVIA
MEDICAL
COMPANY

THE 'PRO'FILE

Clinical Value – The Surgeon's Choice



TIME

A delicate Single-use Instrument that is precise every time

With prompt delivery guaranteed, we endeavour to meet your service needs. Avoiding costly delays, reducing hospital stock holding and improving lead time for Instruments



LIFE

Eliminates risk of contamination from occluded surfaces



COST

Single-use Crocodile Fine Jaw Micro Forceps avoids costly repairs to damaged or worn jaws

This is a cost-effective alternative in comparison to poor quality, cheap Instruments that require opening multiple Packs per patient

No burden of Sterilisation costs

Clinicians' Experiences

Don't just take our word for it ...

"Really happy with the Fine Jaw Crocodile Micro Forceps for Outpatients."

ENT Department, NHS Greater Glasgow and Clyde

See our '9 Things that Make Perfect Single-use Crocodile Forceps' here:
www.dtrmedical.com/products/crocodile-fine-jaw-micro-forceps/

YOUR SINGLE-USE SOLUTION WITH CONFIDENCE

t +44 (0)1792 797910 e marketing@dtrmedical.com w www.dtrmedical.com

