



DTR Medical®

AN
INNOVIA
MEDICAL
COMPANY

THE 'PRO'FILE




0.3mm tip

Successful Dilation with our Double Ended Punctum Dilator



0.4mm tip



140mm

The Design Checklist

- ✓ Our 140mm Double Ended Punctum Dilator is manufactured in the UK to ensure the instrument provides consistent precision every time
- ✓ Used in Punctoplasty to dilate the Lacrimal Punctum for syringing and flushing
- ✓ Designed with high quality polished tips to avoid trauma during entry to the punctum
- ✓ Knurled handle supports improved grip which allows increased accuracy in dilation procedures
- ✓ Well balanced and double ended to suit your procedural needs
- ✓ Clearly labelled 17 x Ø0.4mm and 21 x Ø0.3mm tapered ends to ensure a wide range of dilation
- ✓ Lacrimal Probe and Punctum Seeker also available

Clinicians' Experiences

Don't just take our word for it ...

"Thought our new Dilator was excellent and is now looking to share thoughts with colleagues."

Consultant Ophthalmic Surgeon at Royal Victoria Infirmary

"Delighted with the design and have given them the thumbs up."

St John's Hospital, Lothian

Clinical Value - Quality Instrument to Meet Growing Demand



TIME

With prompt delivery guaranteed, we endeavour to meet your service needs. Avoiding costly delays, reducing hospital stock holding and improving lead time for instruments



LIFE

With the proportion of those aged 65 and older expected to make up ¼ of the UK population by 2045¹, our Punctum Dilator helps meet a growing need with those over 40 most at risk from eye diseases²



COST

Transparency - one price for all
No burden of sterilisation costs
No repair costs for blunt or damaged instruments

1. <https://www.ons.gov.uk/peoplepopulationandcommunity/populationandmigration/populationestimates/articles/overviewoftheukpopulation/mar2017>
2. https://healthcare.utah.edu/healthfeed/postings/2014/03/030514_moran_eye_dilation.php

YOUR SINGLE - USE SOLUTION WITH CONFIDENCE

t +44 (0)1792 797910 e marketing@dtrmedical.com w www.dtrmedical.com

