


CORONET

DMEK EndoGlide™

Evolution of DMEK surgery technique...




PATENT PROTECTED

CORONET

Evolution of DMEK Surgery Technique...

The CORONET™ DMEK EndoGlide™ has been specifically designed to provide controlled graft unfolding with minimal cell loss.



CORONET™ DMEK EndoGlide™ 51-826

KEY FEATURES

- Minimal 2.65mm incision required
- Controlled graft unfolding
- Closed cartridge system protects graft
- Ideal for challenging cases
- Uses controlled pull through paracentesis approach
- Supplied sterile, single use, 5 year shelf life, one per box
- Manufactured in the UK

For best results use with:

The EndoGlide™ Forceps and
EndoGlide™ Support Platform.

Find out more at:

www.networkmedical.co.uk/ophthalmic