

CORONET

DMEK Guarded Trephine Punch

More than just an excellent cut...

Now available
in 10.00mm



PATENT PROTECTED

 **NETWORK**
MEDICAL PRODUCTS

AN
INNOVIA
MEDICAL
COMPANY

CORONET

The CORONET™ DMEK Guarded Trephine Punch offers more than just an excellent cut.

Designed to meet specific demands of DMEK surgery the punch supports instrument position during tissue manipulation by providing an ergonomic and anatomically shaped platform.



CORONET™ DMEK Guarded Trephine Punch 51-855

KEY FEATURES

- Sharp, clean cutting guarded trephine with exclusive CORONET™ Cathedral blade technology
- Pre-set depth to score through Descemet's membrane with a simple press
- Patent protected visualisation windows
- Stable platform base
- Flexible non-restrictive vacuum tubing
- Biocompatible and non-glare blade
- Supplied sterile, single use, 5 year shelf life
- Supplied one per box
- Manufactured in the UK

TREPHINE SIZE (mm)

Size	Code
7.50	51-855-7.50
8.00	51-855-8.00
8.50	51-855-8.50
9.00	51-855-9.00
9.50	51-855-9.50
10.00	51-855-10.00