


NETCELL

Series 5000™ PVA Nasal Pack

The NETCELL™ Series 5000™ PVA Nasal Pack has been clinically proven to reduce trauma and improve patient comfort¹

¹ Source International Rhinology 47,4144. 2009 A randomised controlled trial comparing Rapid Rhino Mannheim and NETCELL Series 5000 packs following routine surgery



Box plots of change in VAS pain score for each pack type. The p value for the difference between these 2 data sets is 0.008.

- Unique patent protected nonstick and soft outer surface prevents tissue ingrowth
- Specifically designed for use after septal, turbinate and rhinoplasty surgery
- Provides gentle, equal pressure within the nasal cavity when stenting postoperative bleeding
- Remains hydrated for up to 7 days without re hydration
- Double compressed for easy insertion (50% slimmer than competitor products)
- Outer coating can be used in conjunction with water or oil based medication
- Supplied sterile, single use only, declared 5 year shelf life

Product Description	Dimension LxTxH	Pack Size	Code
Single hydration nasal packs 6.0cm	6.0 x 1.5 x 2.0cm	Box of 10	15-105
Single hydration nasal pack 8.0cm	8.0 x 1.5 x 2.0cm	Box of 10	15-110
Single hydration nasal pack slimpack 8.0cm	6.0 x 1.0 x 2.0cm	Box of 10	15-111