

Suction Handle-Tubing Recommendations

RECOMMENDED TUBING FOR SUCTION HANDLES				
Tubing (non-sterile)	Tubing Diameters	Spigot	Sterile DSTS Code	Suitable for use Suction Handles
SM405	4.8mm ID x 7.9mm OD	4-6mm	DSTS40506	Aspirating Dissectors Bronchoscopy Frazier Magill Microlaryngeal Oesophageal Poole Sinus Suction Regulator
SM445	6.4mm ID x 9.5mm OD	7-11mm	DSTS44506	All of the above (exc Sinus) plus Barron, Schuknecht and Tapered

Other diameters and lengths available ▶▶▶ please contact our Customer Service team for more information

STERILE TUBING AND SPIGOT STARTER PACKS		
Code	Description	Pack Size
DSTS40506	Silicone Tube 4.8mm ID x 7.9mm OD x 600mm long with Spigot, double packed	10 units
DSTS44506	Silicone Tube 6.4mm ID x 9.5mm OD x 600mm long with Spigot, double packed	10 units



DTR Medical®

DTR MEDICAL SUCTION HANDLE RANGE

Aspirating Dissectors Negus, Freer Hockey Stick & Tonsil

Bladed with suction tip



Barron & Bellucci

3 to 7fg and mini 98mm



Bronchoscopy & Oesophageal

450mm or 550mm 12fg



Frazier

6fg to 12fg, inc. fenestrated and cross fenestrated tip. Available with & without suction control



Glass

Metal variant of traditional instrument



House

Suction & irrigation



Magill

210mm with atraumatic tip



Microlaryngeal

280mm to 330mm, 6 and 9fg



Micro Suction Tubes (Verhoeven) & Regulator

110mm & 120mm 14 to 27g, colour coded



Poole

Standard, neonatal & extra long



Schucknecht

18 to 24swg



Sinus

65 & 90 degree, 9fg & 12 fg with atraumatic tip



Tapered

5, 6 & 7 size with tear drop control



Clearway

Clearway Suction Handle



Yankauer

Metal with 3 & 6 hole rose options



Zoellner & Fine Ends

Combined or separate packs, 18 to 26swg fine ends with metal taper

