

EarPopper®

Middle Ear Inflation Device



The EarPopper® balances pressure in the middle ear by delivering a safe, constant stream of air into the nasal cavity. The EarPopper's constant stream of air opens the Eustachian Tube during the moment of swallowing, clearing the middle ear, relieving negative ear pressure and allowing any fluids to drain. The EarPopper is suitable for adults and children 4+ (with adult supervision).

Treats common middle ear pressure problems without medication or surgery like:

- Eustachian tube dysfunction
- Otitis media with effusion (glue ear)
- Ear fullness caused by air travel or scuba diving

USA Ordering Information

Prescription required

Code	Includes	Quantity/Box
EP-2100	4 AAA Batteries	1
EP-1010	Removable, Single-patient, Disposable Nose Pieces	50

Rest of World Ordering Information

Code	Includes	Quantity/Pack
EP-2100A	4 AAA Batteries	1

Clinically proven to relieve negative ear pressure.



Flying



Fluid in Ear



Diving & Swimming